


SAVANT™ POSTERIOR LUMBAR INTERBODY FUSION SYSTEM

Surgical Technique Guide



CURITEVA®

Table of Contents

SAVANT™ POSTERIOR LUMBAR INTERBODY FUSION SYSTEM


Features and Benefits - PEEK	1
Implant Guide - PEEK	2
Features and Benefits - Titanium Plasma Coated PEEK	3
Implant Guide - Titanium Plasma Coated PEEK	4
Features and Benefits - Titanium	5
Implant Guide - Titanium	6
Instrument Guide	7
Posterior Lumbar Disc Prep Instrument Guide	8
Surgical Technique	10
Product Information.....	15

Disclaimer

The surgical technique shown is for illustrative purposes only. Proper surgical procedure is the responsibility of the medical professional. Please reference the package insert for additional information and system instructions.

SAVANT™ PEEK Posterior Lumbar Interbody Fusion System

The Savant PEEK Posterior Lumbar Interbody Fusion System (PLIF) self-distracting, bulleted tip is offered in a variety of lengths, heights and lordotic angles providing an ideal fit for every patient. The superior and inferior surfaces of the construct have a pattern of teeth designed to provide increased stability and to help prevent movement of the device. The Savant Interbody Fusion System features instruments that are simple and intuitive to use.



PEEK Interbody

Part Number: C304-XXXXYY-Z

XXXX - Interbody Footprint, YY - Interbody Height, Z - Lordosis

- PEEK construction
- Footprints (mm): 22L x 10W, 26L x 10W and 30L x 10W
- Heights: 6 – 16mm, offered in 1mm increments
- Lordosis: 0° and 8°



Savant Features

- Contoured, distracting nose to facilitate placement
- Central and lateral graft windows for biologic seeding
- Multi-directional teeth resist migration
- Strategically positioned tantalum markers for great visualization and placement
- Comprehensive discetomy instruments
- Full range of Trials for proper sizing of implant height
- Convex design offers a more accurate fit to the patient's anatomy
- Robust Inserter engagement features to facilitate implant insertion

SAVANT™ PEEK Posterior Lumbar Interbody Fusion System

PEEK PLIF Interbodies

22mm Length, 10mm Width, 0°	
Height	Catalog Number
6mm	C304-221006-0
7mm	C304-221007-0
8mm	C304-221008-0
9mm	C304-221009-0
10mm	C304-221010-0
11mm	C304-221011-0
12mm	C304-221012-0
13mm	C304-221013-0
14mm	C304-221014-0
15mm	C304-221015-0
16mm	C304-221016-0

26mm Length, 10mm Width, 0°	
Height	Catalog Number
6mm	C304-261006-0
7mm	C304-261007-0
8mm	C304-261008-0
9mm	C304-261009-0
10mm	C304-261010-0
11mm	C304-261011-0
12mm	C304-261012-0
13mm	C304-261013-0
14mm	C304-261014-0
15mm	C304-261015-0
16mm	C304-261016-0

30mm Length, 10mm Width, 0°	
Height	Catalog Number
6mm	C304-301006-0
7mm	C304-301007-0
8mm	C304-301008-0
9mm	C304-301009-0
10mm	C304-301010-0
11mm	C304-301011-0
12mm	C304-301012-0
13mm	C304-301013-0
14mm	C304-301014-0
15mm	C304-301015-0
16mm	C304-301016-0


22mm Length, 10mm Width, 8°	
Height	Catalog Number
7mm	C304-221007-8
8mm	C304-221008-8
9mm	C304-221009-8
10mm	C304-221010-8
11mm	C304-221011-8
12mm	C304-221012-8
13mm	C304-221013-8
14mm	C304-221014-8
15mm	C304-221015-8
16mm	C304-221016-8

26mm Length, 10mm Width, 8°	
Height	Catalog Number
8mm	C304-261008-8
9mm	C304-261009-8
10mm	C304-261010-8
11mm	C304-261011-8
12mm	C304-261012-8
13mm	C304-261013-8
14mm	C304-261014-8
15mm	C304-261015-8
16mm	C304-261016-8

30mm Length, 10mm Width, 8°	
Height	Catalog Number
9mm	C304-301009-8
10mm	C304-301010-8
11mm	C304-301011-8
12mm	C304-301012-8
13mm	C304-301013-8
14mm	C304-301014-8
15mm	C304-301015-8
16mm	C304-301016-8

Radiographic Markers

- 3 Pins per interbody – 2 posterior & 1 anterior
- Pin length: 2.5mm; 1mm diameter
- Posterior pins offset .64mm (tooth depth) from endplate (identified in red below)
- Anterior pin centered on midline
- Pins are 1mm from anterior and posterior edges of interbody
- Posterior pins are 2mm from lateral edges of interbody




Lateral Windows

- 6mm x 0° & 7mm x 8° - No Window
- 7 - 11mm - Standard Windows
- 12 - 16mm - Split Windows

SAVANT™ PEEK Posterior Lumbar Interbody Fusion System Featuring Titanium Plasma Coating

The Savant Titanium Plasma Coated Posterior Lumbar Interbody Fusion System (PLIF) self-distracting, bulleted tip is offered in a variety of lengths, heights and lordotic angles providing an ideal fit for every patient. The superior and inferior surfaces of the construct have a pattern of teeth designed to provide increased stability and to help prevent movement of the device. The Savant Interbody Fusion System features instruments that are simple and intuitive to use.




PEEK Interbody with Coating

Part Number: Part Number: C324-XXXXYY-ZCT

XXXX - Interbody Footprint, YY - Interbody Height, Z - Lordosis

- PEEK construction with Titanium Plasma Coating
- Footprints (mm): 22L x 10W, 26L x 10W and 30L x 10W
- Heights: 6 – 16mm, offered in 1mm increments
- Lordosis: 0° or 8°



Savant Features

- Surface topography produces superior implant contact and construct stability
- Modulus of elasticity between cortical and cancellous bone for optimal load sharing
- Contoured, distracting nose to facilitate placement
- Designed to ensure a high degree of compressive strength and dimensional stability
- Central and lateral graft windows for biologic seeding
- Multi-directional teeth resist migration
- Tantalum radiopaque markers to optimize visibility and placement
- Comprehensive discectomy instruments
- Optimal Radiolucency

SAVANT™ PEEK Posterior Lumbar Interbody Fusion System Featuring Titanium Plasma Coating

Titanium Plasma Coated PEEK PLIF Interbodies

22mm Length, 10mm Width, 0°

Height	Catalog Number
6mm	C324-221006-OCT
7mm	C324-221007-OCT
8mm	C324-221008-OCT
9mm	C324-221009-OCT
10mm	C324-221010-OCT
11mm	C324-221011-OCT
12mm	C324-221012-OCT
13mm	C324-221013-OCT
14mm	C324-221014-OCT
15mm	C324-221015-OCT
16mm	C324-221016-OCT

26mm Length, 10mm Width, 0°

Height	Catalog Number
6mm	C324-261006-OCT
7mm	C324-261007-OCT
8mm	C324-261008-OCT
9mm	C324-261009-OCT
10mm	C324-261010-OCT
11mm	C324-261011-OCT
12mm	C324-261012-OCT
13mm	C324-261013-OCT
14mm	C324-261014-OCT
15mm	C324-261015-OCT
16mm	C324-261016-OCT

30mm Length, 10mm Width, 0°

Height	Catalog Number
6mm	C324-301006-OCT
7mm	C324-301007-OCT
8mm	C324-301008-OCT
9mm	C324-301009-OCT
10mm	C324-301010-OCT
11mm	C324-301011-OCT
12mm	C324-301012-OCT
13mm	C324-301013-OCT
14mm	C324-301014-OCT
15mm	C324-301015-OCT
16mm	C324-301016-OCT

22mm Length, 10mm Width, 8°

Height	Catalog Number
7mm	C324-221007-8CT
8mm	C324-221008-8CT
9mm	C324-221009-8CT
10mm	C324-221010-8CT
11mm	C324-221011-8CT
12mm	C324-221012-8CT
13mm	C324-221013-8CT
14mm	C324-221014-8CT
15mm	C324-221015-8CT
16mm	C324-221016-8CT

26mm Length, 10mm Width, 8°


Height	Catalog Number
8mm	C324-261008-8CT
9mm	C324-261009-8CT
10mm	C324-261010-8CT
11mm	C324-261011-8CT
12mm	C324-261012-8CT
13mm	C324-261013-8CT
14mm	C324-261014-8CT
15mm	C324-261015-8CT
16mm	C324-261016-8CT

30mm Length, 10mm Width, 8°

Height	Catalog Number
9mm	C324-301009-8CT
10mm	C324-301010-8CT
11mm	C324-301011-8CT
12mm	C324-301012-8CT
13mm	C324-301013-8CT
14mm	C324-301014-8CT
15mm	C324-301015-8CT
16mm	C324-301016-8CT

Radiographic Markers

- 3 Pins per interbody – 2 posterior & 1 anterior
- Pin length: 2.5mm; 1mm diameter
- Posterior pins offset .64mm (tooth depth) from endplate (identified in red below)
- Anterior pin centered on midline
- Pins are 1mm from anterior and posterior edges of interbody
- Posterior pins are 2mm from lateral edges of interbody




Lateral Windows

- 6mm x 0° & 7mm x 8° - No Window
- 7 - 11mm - Standard Windows
- 12 - 16mm - Split Windows

SAVANT™ Titanium Posterior Lumbar Interbody Fusion System Featuring TRU-LOK® Technology

The Savant Titanium Posterior Lumbar Interbody Fusion System (PLIF) features TRU-LOK®. The TRU-LOK® surface technology results in a proprietary micro-textured surface intended to provide an improved osteogenic environment. The micro-scale roughness of the TRU-LOK® surface technology is intended to reduce the potential for implant migration and to encourage enhanced osseointegration. The interbody is offered in a variety of lengths, heights and lordotic angles offering an ideal fit for every patient.



Titanium Interbody with TRU-LOK®

Part Number: Part Number: C314-XXXXYY-Z

XXXX - Interbody Footprint, YY - Interbody Height, Z - Lordosis

- Titanium construction with surface texturing
- Footprints (mm): 22L x 10W, 26L x 10W and 30L x 10W
- 6 – 16mm, offered in 1mm increments
- Lordosis: 0° and 8°



Savant Features

- Contoured, distracting nose to facilitate placement
- Oversized central and lateral graft windows for biologic seeding
- Micro-roughness surface texturing produces superior implant contact and construct stability
- Multi-directional teeth resist migration
- Designed to ensure a high degree of compressive strength and dimensional stability
- Open Architecture
- Comprehensive discectomy instruments
- Full range of Trials for proper sizing of implant height

SAVANT™ Titanium Posterior Lumbar Interbody Fusion System Featuring TRU-LOK® Technology

Titanium PLIF Interbodies with TRU-LOK Technology

22mm Length, 10mm Width, 0°	
Height	Catalog Number
6mm	C314-221006-0
7mm	C314-221007-0
8mm	C314-221008-0
9mm	C314-221009-0
10mm	C314-221010-0
11mm	C314-221011-0
12mm	C314-221012-0
13mm	C314-221013-0
14mm	C314-221014-0
15mm	C314-221015-0
16mm	C314-221016-0

26mm Length, 10mm Width, 0°	
Height	Catalog Number
6mm	C314-261006-0
7mm	C314-261007-0
8mm	C314-261008-0
9mm	C314-261009-0
10mm	C314-261010-0
11mm	C314-261011-0
12mm	C314-261012-0
13mm	C314-261013-0
14mm	C314-261014-0
15mm	C314-261015-0
16mm	C314-261016-0

30mm Length, 10mm Width, 0°	
Height	Catalog Number
6mm	C314-301006-0
7mm	C314-301007-0
8mm	C314-301008-0
9mm	C314-301009-0
10mm	C314-301010-0
11mm	C314-301011-0
12mm	C314-301012-0
13mm	C314-301013-0
14mm	C314-301014-0
15mm	C314-301015-0
16mm	C314-301016-0

22mm Length, 10mm Width, 8°	
Height	Catalog Number
7mm	C314-221007-8
8mm	C314-221008-8
9mm	C314-221009-8
10mm	C314-221010-8
11mm	C314-221011-8
12mm	C314-221012-8
13mm	C314-221013-8
14mm	C314-221014-8
15mm	C314-221015-8
16mm	C314-221016-8

26mm Length, 10mm Width, 8°	
Height	Catalog Number
8mm	C314-261008-8
9mm	C314-261009-8
10mm	C314-261010-8
11mm	C314-261011-8
12mm	C314-261012-8
13mm	C314-261013-8
14mm	C314-261014-8
15mm	C314-261015-8
16mm	C314-261016-8

30mm Length, 10mm Width, 8°	
Height	Catalog Number
9mm	C314-301009-8
10mm	C314-301010-8
11mm	C314-301011-8
12mm	C314-301012-8
13mm	C314-301013-8
14mm	C314-301014-8
15mm	C314-301015-8
16mm	C314-301016-8

Lateral Windows

- 6 - 7mm x 0° - No Window
- 7 - 8mm x 8° - No Window
- 8 - 16mm x 0° - Single Window
- 9 - 16mm x 8° - Single Window

Instrument Guide




C305-100

TLIF/PLIF Inserter




C305-500

Straight Tamp



C305-700 PLIF Graft Packing Block



C305-1000 T-Handle, 1/4", Non-Ratcheting



C305-300	6mm x 10mm PLIF Parallel Trial
C305-301	7mm x 10mm PLIF Parallel Trial
C305-302	8mm x 10mm PLIF Parallel Trial
C305-303	9mm x 10mm PLIF Parallel Trial
C305-304	10mm x 10mm PLIF Parallel Trial
C305-305	11mm x 10mm PLIF Parallel Trial
C305-306	12mm x 10mm PLIF Parallel Trial
C305-307	13mm x 10mm PLIF Parallel Trial
C305-308	14mm x 10mm PLIF Parallel Trial
C305-309	15mm x 10mm PLIF Parallel Trial
C305-310	16mm x 10mm PLIF Parallel Trial



C305-350	7mm x 10mm PLIF Lordotic Trial
C305-351	8mm x 10mm PLIF Lordotic Trial
C305-352	9mm x 10mm PLIF Lordotic Trial
C305-353	10mm x 10mm PLIF Lordotic Trial
C305-354	11mm x 10mm PLIF Lordotic Trial
C305-355	12mm x 10mm PLIF Lordotic Trial
C305-356	13mm x 10mm PLIF Lordotic Trial
C305-357	14mm x 10mm PLIF Lordotic Trial
C305-358	15mm x 10mm PLIF Lordotic Trial
C305-359	16mm x 10mm PLIF Lordotic Trial

Posterior Lumbar Disc Prep Instrument Guide



C445-100 8mm Straight Rasp, Bayoneted



C445-125 4mm Straight Curette, Bayoneted



C445-105 8mm Angled Left Rasp, Bayoneted



C445-130 4mm Angled Left Curette, Bayoneted



C445-110 8mm Angled Right Rasp, Bayoneted



C445-135 4mm Angled Right Curette, Bayoneted



C445-115 10mm Straight Osteotome, Bayoneted




C445-140 4mm Push Curette, Bayoneted



C445-120 10mm Curved Osteotome, Bayoneted



C445-145 4mm Pull Curette, Bayoneted



C445-210 Slotted Mallet



C445-150 7mm Ring Curette, Bayoneted



C445-155 7mm Rake Curette, Bayoneted



C445-220 Slap Hammer



C445-160 12mm Cobb Elevator, Bayoneted



C445-165 Nerve Root Retractor, Bayoneted

Posterior Lumbar Disc Prep Instrument Guide



C445-800

Shaver / Distractor, 6mm



C445-801

Shaver / Distractor, 7mm



C445-802

Shaver / Distractor, 8mm



C445-803

Shaver / Distractor, 9mm



C445-804

Shaver / Distractor, 10mm



C445-805

Shaver / Distractor, 11mm



C445-806

Shaver / Distractor, 12mm



C445-807

Shaver / Distractor, 13mm



C445-808

Shaver / Distractor, 14mm



C445-809

Shaver / Distractor, 15mm



C445-810

Shaver / Distractor, 16mm



SAVANTTM

Surgical Technique

Disclaimer

Each surgical step applies to all material types but the illustrations only show PEEK.

Step 1: Patient Positioning

Place the patient in a prone position on a radiolucent operating table that promotes suitable exposure and re-stores sagittal alignment. A/P and lateral fluoroscopy should be used to identify and target the appropriate levels. **(Figure 1)**




Figure 1

Step 2: Surgical Approach to the Disc

A midline incision provides exposure of the interlaminar space and facet joints at the indicated level. Pedicle screws of the surgeon's preference are inserted into the pedicles of the vertebrae adjacent to the disc space to be fused. Using a combination of surgical instruments appropriately selected by the surgeon, a laminotomy, laminectomy and/or facetectomy, is performed, along with the removal of the ligamentum flavum, to gain access to the disc space and identify neural and bony anatomy. **(Figure 2)**




Figure 2

Step 3: Discectomy and Endplate Preparation

Perform a standard discectomy using Paddle Shavers, Osteotomes, Curettes or Rasps to prepare the implant bed. (Figure 3a, 3b)

Note: Adequate preparation of the endplates is important to facilitate vascularization of the bone graft.




Figure 3a




Figure 3b

Step 4: Distraction

Proper distraction is essential to restore desired disc height. Assemble the desired size of the Paddle Shaver/Distractor onto the T-Handle and insert into the disc space. (Figure 4a) Rotating the Paddle Shaver/Distractor clockwise will shave the endplates (Figure 4b) and rotating the Paddle Shaver/Distractor counterclockwise 90° will provide blunt distraction. (Figure 4a)

Note: The Distractors are depth marked from 20 – 40mm to aid visualization.

Note: There are length indicators visible under fluoroscopy

Figure 4a





Figure 4b



Step 5: Implant Selection

Assemble the desired Trial onto the T-Handle. Trials are 10mm wide and 30mm long with length indicators that are visible with fluoroscopy. Heights start at 6mm and increase sequentially by 1mm increments. If necessary, impact the Trial assembly to confirm the position and fit of the Trial. Use A/P and lateral fluoroscopy to confirm proper placement and trajectory. Once the proper size is determined, remove the Trial from the intervertebral space.

A Slap Hammer is available in the Disc Preparation Set for assisting in removing the Trial assembly from the Disc Space. **(Figure 5a, 5b)**

Figure 5a





Figure 5b



Step 6: Implant Preparation

Select the appropriate implant size and attach to the Inserters. Load implant onto the inserter by aligning anti-rotation features and threading central stylus into implant. Once the implant is fully engaged on the Inserter, fill the implant with autologous bone, allograft or other bone grafting material. Repeat this step for the second implant. **(Figure 6a, 6b)**

Alternatively, utilize the Graft Block and Tamp to pack the implant with grafting material. **(Figure 6c)**

Note: The graft should be flush with the upper and lower surfaces of the implant.

Figure 6a




Figure 6b





Figure 6c




Step 7: Implant Insertion

Prior to implant insertion, bone graft material may be placed anteriorly or contralaterally. Under fluoroscopy, insert the implant mounted on the Inserter into the disc space. Once in position, remove all instruments and fill the remaining disc space with bone graft.

Note: A Tamp is available for tamping the implant into final position. (Figure 7b)

To disengage the implant from the Inserter, unthread the knob on the end of the Inserter. (Figure 7a, 7b, 7c)



Repeat steps 3-7 for contralateral implants.

Step 8: Implant Removal (optional)

Attach inserter into implant posteriorly and thread inserter into threads of implant until tightened. Gently remove implant from disc space. If the implant cannot be easily removed, a Cobb elevator, Slap Hammer or Osteotome should be used to loosen the bone to implant interface. (Figure 8)




Figure 8

Indications for Use:

The Savant Lumbar Interbody Fusion System is indicated for use in skeletally mature patients with degenerative disc disease (DDD) at one or two contiguous levels from L2 – S1. DDD is defined as back pain of discogenic origin with degeneration of the disc confirmed by patient history and radiographic studies. These DDD patients may also have up to Grade I spondylolisthesis or retrolisthesis at the involved level(s). Implants are intended to be used with autograft and/or allograft bone (comprised of cancellous and/or corticocancellous bone graft) and supplemental spinal fixation systems that have been cleared for use in the lumbar spine. Patients should receive at least six (6) months of non-operative treatment prior to treatment with the device.

Contraindications:

Contraindications for the Savant Lumbar Interbody Fusion System are comparable to those of other systems of similar design, and include, but are not limited to:

- Patients with probable intolerance to the materials used in the manufacture of this device.
- Patients with infection, inflammation, fever, tumors, elevated white blood count, obesity, pregnancy, mental illness and other medical conditions which would prohibit beneficial surgical outcome.
- Patients unwilling or unable to follow post-operative restrictions on movement, especially in athletic and occupational activities.
- Use with components from other systems, or in any case requiring the mixing of metals from different components.
- Grossly distorted anatomy caused by congenital abnormalities.
- Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery.
- Rapid joint disease, bone absorption, osteopenia. Osteoporosis is a relative contraindication since this condition may limit the degree of obtainable correction, stabilization, the amount of mechanical fixation, and/or the quality of the bone graft.
- Any case where the implant components selected for use would be too large or too small to achieve a successful result.
- Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality.
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance.
- Any case not described in the indications for use.
- Reuse or multiple uses.

Cautions, Precautions and Warnings:

Cautions:

Mixing of dissimilar metals can accelerate the corrosion process. Stainless steel and titanium components must NOT be used together.

Do not use components of the Savant Lumbar Interbody Fusion System with components from any other manufacturer.

Care must be taken to protect the components from being marred, nicked or notched as a result of contact with other objects. Alterations will produce defects in surface finish and internal stresses which may become the focal point for eventual breakage of the implant.

As with all orthopedic implants, none of the Savant Lumbar Interbody Fusion System components should ever be reused under any circumstances.

Precautions:

The implantation of properly selected and placed system implants and components should be performed only by experienced spinal surgeons with specific training in the use of this spinal system because this is a technically demanding procedure presenting a risk of serious injury to the patient.

Patients who smoke have been shown to have an increased incidence of non-union. These patients should be advised of this fact and warned of the consequences. Other poor candidates for spine fusion include obese, malnourished, those with poor muscle and bone quality, and nerve paralysis patients.

Due to the presence of implants, interference with roentgenographic, CT and/or MR imaging may result. The Savant Lumbar Interbody Fusion System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Savant Lumbar Interbody Fusion System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Warnings:

This device system is not intended to be the sole means of spinal support. Its use without a bone graft or in cases that develop into a nonunion will not be successful. No spinal implant can withstand the loads of the body without maturation of a solid fusion mass, and in this case, bending, loosening or fracture of the implant will eventually occur. The proper selection and compliance of the patient will greatly affect the results.

The implantation of spinal systems should be performed only by spinal surgeons fully experienced in the surgical techniques required for the use of such implants. Even with the use of spinal implants, a successful result in terms of pain, function, or fusion is not always achieved in every surgical case.

The physician is the learned intermediary between the company and the patient. The indications, contraindications, warnings, and precautions given in this document must be conveyed to the patient. If requested, additional information, including surgical technique manuals, may be obtained through Curiteva customer support representatives.

Toll Free: 877.9CURITEVA
Phone: 256.213.1057
Fax: 256.213.1058
www.curiteva.com
customersupport@curiteva.com



CURITEVA®
POWER TO HEAL

Products patented and/or patent pending.
All products are not currently available in all markets.

Copyright©2020 by Curiteva, Inc.
QMSF-1145 Rev 1.0