



CD HORIZON® SOLERA® Spinal System

Surgical Technique



CAPSTONE® PEEK Spinal System, CLYDESDALE® Spinal System, CRESCENT® Spinal System—PEEK, CRESCENT® Spinal System—Titanium, and SOVEREIGN® Spinal System incorporate technology developed by Gary K. Michelson, MD.



Realign, Reconstruct, Customize.

From the thoracic spine to the ilium, the CD HORIZON® SOLERA® Spinal System facilitates surgeon choice and flexibility across patient types with a variety of implant options for treating multiple spinal pathologies with one system family. Spinal rod size options include 4.75mm, 5.5mm, and 6.0mm diameters; and material options of commercially pure titanium (not available in 4.75mm rod), titanium alloy, CHROMALOY™ and CHROMALOY™ Plus. Risks associated with these spinal implants include loosening, disassembly, bending and/or breakage of components.

The system offers the opportunity to reduce overall metal mass and profile, and provide surgeon choice without compromising implant integrity, as seen in mechanical testing when compared to the CD HORIZON® LEGACY™ System.¹ The technology platform offered with this system is backed by more than 25 years of CD HORIZON® clinical experience and Medtronic expertise.²




Low Profile Implants



Rod Material Choices



Power, Neuromonitoring, and Navigation Compatibility



Minimally Invasive Approach Options



Advanced Instrumentation

Implant Traceability



¹ Based on mechanical testing of the CHROMALOY™ and CHROMALOY™ Plus rod construct per ASTM F1798.

² Based on internal sales estimates.



CD HORIZON[®] SOLERA[®] Spinal System

Surgical Technique

Implant Features	2	Hook Surgical Technique	
Titanium Hook Implants	6	Surgical Strategy	60
Instrument Overview	7	Hook Site Preparation, Options, and Insertion	61
Instrument Set	9	Hook Site Preparation, Options, and Placement	62
Pedicle Screw Surgical Technique		Decortication	66
Thoracic Facetectomy and Starting Points	22	Rod Contouring	67
Screw Starting Points	23	Rod Insertion	68
Screw Options	24	Deformity Correction	69
Pedicle Preparation	25	Compression and Distraction	71
Enabling Technologies	27	Stabilization and Holding Rod Placement	72
VERIFY [®] Implant Tracking System	28	Final Tightening	73
Screw Placement	29	CD HORIZON [®] X10 CROSSLINK [®]	74
ATS [™] AWL TAP Multi-Axial Screw Implant		Plate Placement and Closure	
Placement Technique	31	Suggestions to Explant the	75
Lateral Connectors	32	CD HORIZON [®] SOLERA [®] Spinal System	
Interbody or Vertebral Body Replacement Options	34	Product Ordering Information	76
Rod Selection	35	Important Product Information	86
Rod Contouring and Placement	38		
Provisional Driver Assembly	39		
Rod Reduction	40		
Sequential Rod Reduction	42		
5.5/6.0mm Lateral Translator	44		
5.5/6.0mm SMARTLINK [™] Instruments	45		
Deformity Correction	50		
Final Tightening and Decortication	53		
Bone Graft Options	54		
CD HORIZON [®] X10 CROSSLINK [®] Plate Placement	55		
Final Tightening of the Plate	57		
Postoperative Care and Mobilization	58		

Implant Features

Multi-Axial Screw (MAS)



- » 4.75mm screw features a 26% reduction in overall volume than CD HORIZON® LEGACY™ 5.5mm MAS and 12% less volume than CD HORIZON® LEGACY™ 4.5mm MAS
- » 5.5/6.0mm screw features a 10% reduction in overall volume than CD HORIZON® LEGACY™ 6.35mm MAS and is compatible with a 5.5mm rod or a 6.0mm rod diameter
- » OSTEOGRIP® dual lead threadform
- » Cobalt Chrome tulip is designed to be compatible with CHROMALOY™, CHROMALOY™ Plus, and Titanium Rods, allowing real-time rod choices in the OR to customize the stiffness and strength of the construct
- » Imaging comparable to Titanium screw heads¹
NOTE: The CD HORIZON® Spinal System has not been evaluated for safety, heating, migration or compatibility in the magnetic resonance environment.
- » Saddle is color-coded by bone screw diameter
- »

ATS™ 4.75mm AWL-TAP Multi-Axial Screw*




Order set number SPS02684.

- » Available for 4.75mm rod diameter systems
- » Awl tip with cutting flute eliminates five procedural steps compared to the common screw placement technique:
 - Less instruments required in the sterile field during screw placement
 - Less instrument exchanges during screw placement
- » Fully threaded dual lead bone screw
- » Use with image guidance is recommended

*Compatible with CD HORIZON® SOLERA® SEXTANT® and CD HORIZON® LONGITUDE® II 4.75mm Extenders.

Sagittal Adjusting Screw (SAS)



- » A fixed pedicle screw that combines the sagittal forgiveness of a poly axial screw and the direct vertebral body control of a mono-axial screw
- » Enables sagittal adjustment of vertebral bodies +/-13° cephalad and caudal
- » Accepts 5.5mm and 6.0mm rod diameters to accommodate construct demand
- » Saddle is color-coded by bone screw diameter

For more information, refer to the CD HORIZON® SOLERA® Spinal System Pathology Surgical Technique.

Screw Color-coding Size Reference

4.0mm	4.5mm	5.0mm	5.5mm	6.0mm	6.5mm	7.5mm	8.5mm	9.5mm	10.5mm

¹Results based on internal imaging study using the CD HORIZON® SOLERA® 4.75mm Spinal System. MRI and CT images were taken of three different materials: Stainless Steel, Titanium, and Cobalt Chrome. Images were reviewed by seven technical experts for clarity in regions of interest.

Implant Features *continued*

Fixed Angle Screw (FAS)



- » 4.75mm FAS is 21% smaller than CD HORIZON® LEGACY™ System 5.5mm Fixed Angle Screw
- » 5.5/6.0mm FAS features a 10% reduction in overall volume than CD HORIZON® LEGACY™ 6.35mm FAS and is compatible with a 5.5mm rod or a 6.0mm rod diameter
- » OSTEOGRIP® dual lead threadform
- » Color-coded by bone screw diameter
- » Compatible with CHROMALLOY™, CHROMALLOY™ Plus, and Titanium Rods, allowing real-time rod choices in the OR to customize the stiffness and strength of the construct

Multi-Axial Reduction Screw (MARS)



- » Available for both 4.75mm and 5.5/6.0mm rod diameter systems
- » Like the Multi-Axial Screw, these screws can be inserted with the bone screw shaft at the most convenient angle of insertion
- » Extended tabs allow for a larger window to capture and slowly reduce the spinal rod and are broken off when rod reduction is complete

For more information, refer to the CD HORIZON® SOLERA® System Multi-Axial Reduction Screw Surgical Technique.


Hydroxyapatite (HA) Coated Multi-Axial Screw



- » Available for both 4.75mm and 5.5/6.0mm rod diameter systems
- » OSTEOGRIP® dual lead threadform
- » Hydroxyapatite (HA) coating
- » Risks associated with these spinal implants include loosening, disassembly, bending, and/or breakage of components

Implant Features *continued*

Break-off Set Screw



Blunt Start Thread Cut

- » Features a blunt start thread, reducing the chances of the set screw starting off-axis to the tulip
- » The thread pattern on the set screw and the geometry of the tulip head forces the set screw to start in "one way."
- » The reverse angle threadform maximizes the surface contact of the set screw threads with the tulip head
- » Internal mechanical testing shows increased performance with an average of 17% reduction in locking torque of the 4.75mm system vs. the CD HORIZON® LEGACY™ 5.5mm system, and 4% reduction of the 5.5/6.0mm system vs. the CD HORIZON® LEGACY™ 6.35mm System. Mechanical testing is not indicative of clinical results.

Rod Spectrum



CHROMALLOY™ Pre-bent Rod



CHROMALLOY™ Rod



CHROMALLOY™ Plus Rod



Commercially Pure Titanium Rod



Titanium Alloy Rod*



CHROMALLOY™ Plus Apex Rod*



CHROMALLOY™ Plus Tapered Rod*

- » Multiple material types and rod diameters for construct tailoring and intraoperative flexibility
- » Available rods
 - Precontoured and Precut CHROMALLOY™ and Commercially Pure Titanium Rods
 - Straight CHROMALLOY™, CHROMALLOY™ Plus, Commercially Pure Titanium, and Titanium Alloy Rods
 - Straight CHROMALLOY™ Plus APEX Rods
 - Straight CHROMALLOY™ Plus Tapered Rods
 - See "Product Ordering Information" on page 79 and 80 for the comprehensive list of rod options

*(can be ordered as an extra implant)

Implant Features *continued*

CD HORIZON® X10 CROSSLINK® Plate



- » Available for 4.75mm, 5.5mm, and 6.35mm rod diameter systems
- » Features a reduced profile to better match the lower profile Multi-Axial Screw
- » Compatible with CHROMALLOY™, CHROMALLOY™ Plus, and Titanium Rods
- » Adjustable or fixed length options
- » Adjustable plate attaches to the rod in the coronal, sagittal, or transverse planes in any orientation

Hooks



- » Compatible with CHROMALLOY™, CHROMALLOY™ Plus, and Titanium Rods
- » Anatomic design to mimic the posterior spinal elements

Closed Multi-Axial Screws (CMAS)




- » Closed Multi-Axial Screws (CMAS) are compatible with the 6.0mm Lateral Connector Post and will provide flexibility to position the screw head



Compatible Set Screw 7049855

Lateral Connectors



- » Compatible with CHROMALLOY™, CHROMALLOY™ Plus, and Titanium Rods
- » Allows medial/lateral adjustability which may minimize the need for coronal rod bending



Compatible Set Screw 779170005



Compatible Break-off Set Screw
4.75mm 5440030
5.5/6.0mm 5540030

Titanium Hook Implants

	Hook Type	Vertebral Posterior Element Placement	Blade Direction	Region of Spine	Design Features
	Pedicle Hook	Articular Process	↑	T1 – T10	» Bifid blade grasps thoracic pedicle for stability.
	Wide Blade Hook	Lamina	↕	T1 – L5	» Wider blade width distributes forces evenly over a wider aspect of bone.
		Transverse Process	↕	T1 – L5	
	Narrow Blade Hook	Lamina	↕	T1 – L5	» Narrower blade width minimizes metal volume in the spinal canal.
		Transverse Process	↕	T1 – L5	
	Wide Blade Ramped Hook	Lamina	↕	T1 – L5	» Ramp reduces intra-canal intrusion.
		Transverse Process	↕	T1 – L5	
	Narrow Blade Ramped Hook	Lamina	↕	T1 – L5	» Ramp reduces intra-canal intrusion.
		Transverse Process	↕	T1 – L5	
	Extended Body Hook	Lamina	↕	T1 – L5	» Can correct anatomic misalignment between two laminae in the dorso-ventral plane.
		Transverse Process	↕	T1 – L5	
	Offset Hook	Lamina	↕	T1 – L5	» Can be used to medialize or lateralize the rod in supralaminar or infralaminar position. » Can back up a pedicle screw at the same level.
		Transverse Process	↕	T1 – L5	
	Total Anatomical Pedicle Hook	Articular Process	↑	T1 – T10	» Centralized head for balance. » Lipped design can improve hook stability.
	Total Anatomical Transverse Process Hook	Transverse Process	↕	T1 – L5	» Centralized head for balance. » Lipped design can improve hook stability.
	Total Anatomical Supralaminar Hook	Lamina	↓	T1 – T10	» Upward slope of blade and body approximates slope of the lamina.
	Total Anatomical Infralaminar Hook	Lamina	↑	T10 – L2	» Downward slope of blade and body approximates the inferior margin of the lamina.

Color-coding Size Reference

Extra Small	
Small	
Medium	
Large	


Improper positioning and placement of implants may cause damage.

Instrument Overview

5.5/6.0mm SMARTLINK™ Extenders



5.5/6.0mm Lateral Translator



4.75mm Sequential Reducer



Dual Ended Set Screw Starter

Securely engages the outer hex of the set screw to allow more control when tightening the set screw.



Instrument Overview *continued*

Hinged Translator




Compresses and distracts off of the rod, not the implant.

Rod Gripper



Attaches firmly to the rod for controlled rod manipulation.

Adjustable Handle



Allows easy attachment to the Provisional Driver and Rod Pusher shafts. The short and extended positions allow for customization of preferred lengths.

Instrument Set


Screw Preparation



In-Line Handle Round Awl
7480104



Lumbar Probe
8350293



Thoracic Probe
8350294



Quick Connect Ratcheting Handle
G900000



Quick Connect Ratcheting T-Handle
7579000



Taps
8350420 (3.75mm)
8350421 (4.0mm)
8350422 (4.5mm)
8350423 (5.0mm)
8350424 (5.5mm)
8350425 (6.0mm)
8350426 (6.5mm)
8350428 (7.5mm)
8350430 (8.5mm)
8350432 (9.5mm)

Dual Lead Taps*
5480035 (3.75mm)
5480040 (4.0mm)
5480045 (4.5mm)
5480050 (5.0mm)
5480055 (5.5mm)
5480060 (6.0mm)
5480065 (6.5mm)
5480075 (7.5mm)
5480085 (8.5mm)




Dual Ended Feeler Probe
7480100



Sounding/Feeler Probe, 2mm
8572102



Straight Holt Probe*
803-292



T-handle Holt Probe*
803-290

*May be ordered as an extra instrument.

Instrument Set *continued*

Screw Placement



Fixed Angle Screw Lock Sleeve Driver

5484334 – 4.75mm
5584334 – 5.5/6.0mm



Lock Sleeve Multi-Axial Screwdriver

5484306 – 4.75mm
5584111 – 5.5/6.0mm



Multi-Axial Screwdriver*

5484305 – 4.75mm
5584109 – 5.5/6.0mm

*May be ordered as an extra instrument.

Screw Placement *continued*



Self-Retaining Driver

5484113



Screw Confirmation Tool


5584007

Hook Preparation



Laminar Elevator

7480220



Pedicule Finder

7480225



Transverse Process Elevator

7480222

Instrument Set *continued*

Hook Placement



Hook Pusher for Implant Holder
5484231 – 4.75mm
5584418 – 5.5/6.0mm



Self-Retaining Implant Holder
5484217 – 4.75mm



Lateral Implant Holder
5484211 – 4.75mm
7480211 – 5.5mm



Offset Hook Holder
5484219 – 4.75mm
5584119 – 5.5/6.0mm

Hook Placement *continued*



Straight Implant Holder
5484216 – 4.75mm
7480216 – 5.5mm*



Self-Retaining Hook Pusher
5484317 – 4.75mm
5584317 – 5.5/6.0mm

*May be ordered as an extra instrument.

Instrument Set *continued*

Rod Contouring



50cm Rod Template
5484510



20cm Rod Template
5484520



4.75mm Rod Cutter
RC404022-E

Table Top Rod Cutter
808-529



French Bender
7480162

Rod Insertion



Rod Introducer/Holder
5484318 – 4.75mm
5584312 – 5.5/6.0mm

Instrument Set *continued*

Rod Reduction

**Beale Rod Reducer, Slots**5484307 – 4.75mm
5584134 – 5.5/6.0mm**Forceps Rocker**

7480142 – 4.75mm

**Sequential Reducer Inner Sleeve**


5484304 – 4.75mm

**Push Button Rocker***

5484319 – 4.75mm

**Sequential Reducer Outer Sleeve**

5484328 – 4.75mm

**4.75mm Sequential Reducer
(Assembled)****Rod Pusher**5484235 – 4.75mm
7486235 – 5.5/6.0mm

Shown with Adjustable Driver Handle

*May be ordered as an extra instrument.

Instrument Set *continued*

Rod Reduction *continued*



Lateral Translator
5485904 – 5.5/6.0mm



Translator Driver
5485905 – 5.5/6.0mm



Open Extender
5485900 – 5.5/6.0mm



Closed Extender
5485901 – 5.5/6.0mm



Extender Guide
5485906 – 5.5/6.0mm



Reducer, HS
5485908 – 5.5/6.0mm




Driver
5485903 – 5.5/6.0mm



SMARTLINK™ Extender
(Assembled)

Instrument Set *continued*

Derotation



Global Link
(5485913)




Interlink
(5485911)



Lateral Handle
(5485914)



Screw Locker*
(5485912)



Segmental Link
(5485910)

*May be ordered as an extra instrument.


Instrument Set *continued*

Rod Correction



Coronal Benders

Left	Right
5484315 – 4.75mm	5484316 – 4.75mm
5584265 – 5.5/6.0mm	5584270 – 5.5/6.0mm



In Situ Benders

Left	Right
5484313 – 4.75mm	5484314 – 4.75mm
5584255 – 5.5/6.0mm	5584260 – 5.5/6.0mm



Dual Action Rod Gripper

5484303 – 4.75mm
5584006 – 5.5/6.0mm



Rod Rotation Wrench

5484285 – 4.75mm
7480285 – 5.5/6.0mm

Instrument Set *continued*

Compression and Distraction



Multilevel Hook Compressor*
808-503



Distractor Pliers
5484171 – 4.75mm
5584171 – 5.5/6.0mm




Curved Compressor
5484309 – 4.75mm
5584165 – 5.5/6.0mm



Right




Left




Small Parallel Compressor
7480165

Hinged Translator
Left 5484311 – 4.75mm 5584173 – 5.5/6.0mm
Right 5484310 – 4.75mm 5584172 – 5.5/6.0mm



Large Parallel Compressor
7480166



TLIF Distractor
5484327 – 4.75mm
5584327 – 5.5/6.0mm



Parallel Distractor
7480170


*May be ordered as an extra instrument.

Instrument Set *continued*

Final Tightening/Set Screws



Dual Ended Set Screw Starter
5484302



Set Screw Breakoff Driver
8350316



Counter Torque
5484308 – 4.75mm
5584150 – 5.5/6.0mm



Torque Indicating Driver
5484147 – 4.75mm
5485147V – 5.5/6.0mm



6.35mm Hex Shaft Provisional Driver
5480131V
Shown with Adjustable Driver Handle




Adjustable Driver Handle
5480140V



Provisional Driver
7480131

Instrument Set *continued*

Removal Instruments



**Ball-ended T25 Bone Screw
Removal Driver**
5480004V




Self-Retaining 3.5mm Hex Driver
7480114




T25 Shaft
815-518



3.0mm Hex Driver
8110530



T27 Obturator
7480154




T25 Set Screw Obturator
8350322

Additional Instruments



Torx 20 Shaft
815-517



**Screw Head
Positioner Shaft**
5484301 – 4.75mm
5584003 – 5.5/6.0mm



Pedicule Marker, Cylindrical
8360912



Pedicule Marker with Hex
8360911



Plate Holding Pin Driver
9790903

Instrument Set *continued*

CD HORIZON® X10 CROSSLINK® Plate Instrument Set




Measuring Card
5484330 – 4.75mm
8110501 – 5.5/6.35mm



T-bolt Implant Positioners
808-545



Measuring Caliper
5484329 – 4.75mm
8110502 – 5.5/6.35mm



7/32" Torque-Limiting Set Screwdriver
8110535



In-line Plate Holder
8110511




Counter Torque
8110540



3.0mm Hex Head Shaft, Removal Driver
8110530



Plate Benders Bending Irons
8114505 – 4.75mm
8110525 – 5.5/6.35mm



Pedicle Screw Surgical Technique



Thoracic Facetectomy and Starting Points

Clean the facet joints and perform a partial inferior articular process osteotomy to enhance visualization and fusion. Remove 3mm to 5mm of the inferior facet and denude the articular cartilage on the superior facet, except for the lowest vertebra to be instrumented. This will allow for the intraoperative localization of the thoracic pedicle screw starting points.

Anatomic starting points vary by the posterior elements that can be observed intraoperatively. These include the transverse process, the lateral portion of the pars interarticularis, and the base of the superior articular process.

After a thorough exposure, use as much anatomic information as possible by starting with a neutral, non-rotated vertebra. The lateral and posterior views shown on the following page can be used as a guide for starting points and screw trajectory.

The first and extremely critical step to performing advanced deformity techniques is the safe and secure placement of segmental pedicle screws. Knowledge of the Superior Facet Rule (A) to direct the medial/lateral and the Cephalocaudal Starting Points (B) is a helpful reference to accomplish this.




! Important

Do not start medial to the midpoint of the superior facet.


Screw Starting Points

Use Fixed Angle or Multi-Axial Screws for the straightforward approach (Blue Pins).
Use Multi-Axial Screws only for the anatomic approach (Green Pins).



Axial View

Level	Cephalad-Caudad Starting Point	Medial-Lateral Starting Point
T1	Midpoint Transverse Process (TP)	Junction: TP-Lamina
T2	Midpoint TP	Junction: TP-Lamina
T3	Midpoint TP	Junction: TP-Lamina
T4	Junction: Proximal Third-Midpoint TP	Junction: TP-Lamina
T5	Proximal Third TP	Junction: TP-Lamina
T6	Junction: Proximal Edge-Proximal Third TP	Junction: TP-Lamina-Facet
T7	Proximal TP	Midpoint Facet
T8	Proximal TP	Midpoint Facet
T9	Proximal TP	Midpoint Facet
T10	Junction: Proximal Edge-Proximal Third TP	Junction: TP-Lamina-Facet
T11	Proximal Third TP	Just medial to lateral pars
T12	Midpoint TP	At the level of lateral pars
L1	Midpoint TP	Junction: Lateral pars and superior facet
L2	Midpoint TP	Junction: Lateral pars and superior facet
L3	Midpoint TP	Junction: Lateral pars and superior facet
L4	Midpoint TP	Junction: Lateral pars and superior facet
L5	Midpoint TP	Junction: Lateral pars and superior facet
S1	Midpoint Sacral Ala	Junction: Sacral ala and superior facet
Iliac	1cm Cephalad to Distal Posterior Superior Iliac Spine (PSIS)	1cm inferior to the superior PSIS on the medial slope



Oblique View

Screw Options

There are two multi-axial bone screw thread form options offered with the CD HORIZON® SOLERA® 4.75mm Spinal System and the preparation and placement techniques differ. The chart below shows a side-by-side comparison of each procedure. The follow pages describe standard screw placement. For ATS™ AWL TAP Multi-Axial screw placement steps, please see page 31.

Traditional Implant Procedure

1. Create a 3mm-deep posterior cortical breach with a high-speed burr.
2. Insert a probe to the appropriate depth and then rotate 180 degrees to make room for the screw.
3. Use a flexible ball-tipped Sounding/Feeler Probe to confirm five distinct bony borders: a floor and four walls (medial, lateral, superior, and inferior).
4. Undertap the pedicle by 0.5mm to 1.0mm of the final screw diameter.
5. Palpate the tapped pedicle tract with a flexible Sounding/Feeler Probe.
6. Clamp a hemostat to the exposed Sounding/Feeler Probe and measure the length of the hole.
7. Select the appropriate screw diameter and length by both preoperative measurement and intraoperative observation.
8. Attach the standard Multi-Axial Screw to the Multi-Axial Screw Lock Sleeve Driver.
9. Slowly advance the screw down the pedicle to ensure proper tracking while allowing for viscoelastic expansion.

ATS™ Implant Procedure

1. Access the pedicle and select the appropriate ATS™ diameter and length by both preoperative measurement and intraoperative observation.
2. Attach the Multi-Axial ATS™ to the Multi-Axial Screw Lock Sleeve Driver.
3. Slowly advance the ATS™ down the pedicle to ensure proper tracking while allowing for viscoelastic expansion.

Pedicle Preparation

Create a 3mm-deep posterior cortical breach with a high-speed burr. A pedicle blush may be visualized suggesting entrance into the cancellous bone at the base of the pedicle. Occasionally, when preparing small pedicles located at the apex of the curve, the blush will not be evident due to the limited intrapedicular cancellous bone. In this case, use the Thoracic Probe to search in the burred cortical breach for the soft, funnel-shaped cancellous bone, which indicates the entrance to the pedicle. The tip should be pointed laterally to avoid perforation of the medial cortex (**Figure 1**).

Grip the side of the handle to avoid applying too much ventral pressure. Insert the tip approximately 20mm to 25mm (**Figure 2**), and then remove the probe to reorient it so that the tip points medially. Carefully place the probe into the base of the prior hole and use the instrument markings to advance the probe to the desired depth (**Figure 3**). Rotate the probe 180° to ensure adequate room for the screw.



Figure 1



Figure 2



Figure 3

Pedicle Preparation *continued*

Check to ensure that only blood is coming out of the pedicle and that the bleeding is not excessive. Using a flexible ball-tipped Sounding/Feeler Probe, advance the instrument to the base (floor) of the hole to confirm five distinct bony borders: a floor and four walls (medial, lateral, superior, and inferior) (**Figure 4**). Give special care to the first 10mm to 15mm of the tract. Cortically breached pedicles may be salvageable. When necessary, place bone wax in the pedicle hole to limit bleeding, then reposition the probe with a more appropriate trajectory.

Next, undertap the pedicle by 0.5mm to 1.0mm of the final screw diameter (**Figure 5**). Palpate the tapped pedicle tract with a flexible Sounding/Feeler Probe. Clamp a hemostat to the exposed Sounding/Feeler Probe and measure the length of the hole (**Figure 6**). Select the appropriate screw diameter and length by both preoperative measurement and intraoperative observation.

NOTE: Dual Lead taps can be ordered as an extra instrument. Although the OSTEGRIP® thread design is different than previous designs, you can still select a tapping strategy based upon personal preference. In most instances we recommend under tapping by 1mm. If a pedicle seems particularly tenuous or brittle, then line-to-line tapping might be considered. In addition, due to the untapped second cortical thread, line to line tapping generates a significantly better sense of rigid fixation with the CD HORIZON® SOLERA® screw as compared to prior systems. For example, you have the flexibility to follow a 5.5mm tap, which feels tight, with a 5.5mm screw and yet not compromise on fixation.



Figure 4



Figure 5



Figure 6


Enabling Technologies

Triggered intraoperative EMG monitoring, such as the NIM-ECLIPSE® Spinal System, may be used to verify the trajectory within the pedicle. The O-ARM® Imaging System coupled with the STEALTHSTATION® Image Guidance System can also be used to navigate pedicle preparation and screw placement.

The POWEREASE® System is compatible with 4.75mm and 5.5/6.0mm CD HORIZON® SOLERA® System implants. The POWEREASE® System includes a rod cutter, post cutter, and set screw break-off instruments that result in reduced physical fatigue for surgeons as compared to manual instruments (based on biomechanical testing and claims validation questionnaire using Likert scale and completed by nine surgeons).



NIM-ECLIPSE® Spinal System



NIM® Pedicle Probes



NAVLOCK®
Navigated Instruments



POWEREASE® System



O-ARM® Imaging System



O-ARM® and STEALTHSTATION® System Images

The NIM-ECLIPSE® Spinal System is manufactured by Medtronic Xomed, Inc. Distributed by Medtronic Spinal. For the complete labeling for the navigation products please contact Medtronic Navigation, General Business at (888) 580-8860 or visit www.medtronicnavigation.com.

VERIFYI® Implant Tracking System

CD HORIZON® SOLERA® Spinal System implants feature the VERIFYI® Implant Tracking System. Each implant has a tag attached that is clearly marked with the part number, lot number, implant size, and a barcode (Figure 7). Contact your local Medtronic sales representative for detailed information on using the VERIFYI® Implant Tracking System.

Prior to implantation remove the tags from the implants (Figure 8). Retain all of the implant tags so that they can be scanned at the end of the surgery. A tag sorter is available if the surgeon wants to track the implants by spinal level (Figure 9).

! Important

The implant tags must be removed prior to implantation. Do not implant the tags.



Figure 7




Figure 8



Figure 9

Screw Placement

Quick Connect Handle

Attach the Quick Connect Handle to the Multi-Axial Screw Lock Sleeve Driver by snapping into place. A slight rotation of the Quick Connect Handle may be required to fully engage with the driver. To remove the Quick Connect Handle from the driver, press the cap on the handle to disengage the driver (Figure 10a and 10b).



Figure 10a

Figure 10b

Screw Engagement

After the Quick Connect Handle is assembled on the Lock Sleeve Driver, ensure that the blue locking cap is not engaged with the screwdriver shaft and then thread the driver shaft into the screw from the screw caddy (Figure 11). Slide the blue locking cap toward the screw to engage it with the driver shaft (Figure 12). An audible "click" will confirm engagement. Break off the VERIFY® Implant Tracking Tag as shown (Figure 13) and place the tag in the Tag Sorter.




Figure 11

Figure 12

Important

The Multi-Axial Screw Lock Sleeve Driver locking cap must be disengaged while threading it into the threads of the bone screw tulip head. The locking mechanism of the driver may be damaged if it is advanced into the tulip head with the locking cap engaged.



Figure 13

Screw Placement *continued*

The ring on the Quick Connect Handle determines the direction the screw will be driven by the Multi-Axial Screw Lock Sleeve Driver. Turn the ring clockwise to drive the screw into the pedicle (**Figure 14**). Turning the ring counterclockwise will allow the driver to remove the screw from the pedicle (**Figure 15**).



Figure 14



Figure 15

Slowly advance the screw down the pedicle to ensure proper tracking while allowing for viscoelastic expansion (**Figure 16**). Once the screw is inserted, push the button on the blue locking cap and slide it back toward the handle to disengage the driver (**Figure 17**). Finally, unthread the Lock Sleeve Multi-Axial Screw Driver from the screw.

! Important

If the button on the blue cap is locked and the driver is difficult to disengage, hold the screwdriver sleeve and twist the screw driver shaft counterclockwise. This will allow the button on the locking cap to release and disengage the screwdriver shaft from the screw.



Figure 16

Screws should be placed at every segment on the correction side and every third or fourth level on the stabilizing side. Insert at least two screws at the proximal and distal ends of the planned construct on the stabilizing side. For some pathologies, such as kyphosis and scoliosis, more screws are placed for greater construct rigidity. Screws should be checked radiographically at this time to ensure intraosseous screw placement.



Figure 17

ATS™ AWL TAP Multi-Axial Screw Implant Placement Technique

The 4.75mm ATS™ bone screw design has a tapered awl/tap tip with cutting flutes which obviates the probing and tapping steps prior to screw placement. Due to the sharp tip design, use of intraoperative imaging is recommended. The implant can be used with the O-arm® Imaging System coupled with the STEALTHSTATION® Image Guidance System and the navigation-compatible IPC® POWEREASE® System. The ATS™ bone screw is compatible with triggered intraoperative EMG monitoring, such as the NIM-ECLIPSE® Spinal System, which may be used to verify the trajectory within the pedicle.

When used with the O-arm® Imaging System coupled with the STEALTHSTATION® Image Guidance System, the screw size is measured intraoperatively. The tip of the screw should be 1cm short of the anterior cortical wall of the vertebral body. Localize the pedicle and assess the patient specific anatomy per the surgeon's standard technique. Use the measurements on the navigated probe to verify the screw length as well as the trajectory. Attach the screw to the navigated driver and dock the implant in the bone (**Figure 18a**). Slowly advance the ATS™ implant down the pedicle to ensure proper tracking while allowing for viscoelastic expansion. If using the IPC® POWEREASE® System for screw placement, please refer to the IPC® POWEREASE® System User Manual.

When placing the ATS™ implant manually, under fluoroscopy, it is important to determine the screw lengths preoperatively. Localize the pedicle and assess the patient specific anatomy per the surgeon's standard technique and use the measurements on the probe to verify the screw length as well as the trajectory. Following intraoperative fluoroscopy, attach the ATS™ implant to the MAS driver (**Figure 18b**). Slowly advance the ATS™ implant down the pedicle to ensure proper tracking while allowing for viscoelastic expansion. EMG monitoring, such as the NIM-ECLIPSE® Spinal System, may be used to verify the trajectory within the pedicle.

✓ Note

In pedicles with predominately thick cortical bone, a smaller diameter screw is recommended.

! Important

It is advised that intraoperative imaging be used to aid screw placement, and the screws should not be used for bicortical fixation.

✓ Note

If particularly hard bone is encountered (for example, dense sclerotic bone), it might be helpful to prepare the pedicle by creating a cortical breach at the pedicle entry point as per surgeons standard technique prior to placing ATS implant.




Figure 18a



Figure 18b

Lateral Connectors

Closed, Top-Loading, and Side-loading Lateral Connectors are color-coded by rod channel diameter and connect to 4.75, 5.5, and 6.0mm diameter rods as shown below (Figure 19a). The titanium break-off set screw (779170005) should be used with Side-loading and Closed Connectors. The Top-loading Connectors are compatible with the set screws 5440030 for 4.75mm and 5540030 for 5.5/6.0mm sizes.

Rod Size Compatibility				
4.75mm	5.5mm	6.0mm	Lateral Connectors	Set Screw
				
				
				
				
				
				

Figure 19a

The smooth posts of the connectors are offered in 20mm to 70mm lengths and are compatible with 4.75, 5.5, 6.0, and 6.35mm CD HORIZON® System Screws (Figures 19b, 19c, and 19d).



Figure 19b



Figure 19c



Figure 19d

Lateral Connectors *continued*

When performing iliac fixation, Closed, Top-Loading, and Side-Loading Lateral Connectors may be used to facilitate construct assembly (**Figure 19e**).





Figure 19e



Closed
Lateral
Connector



Top-Loading
Lateral
Connector



Side-Loading
Lateral
Connector

These connectors are compatible with long CD HORIZON® SOLERA™ Multi-Axial Screws and CD HORIZON® Closed Multi-Axial Iliac Screws to allow intraoperative flexibility (**Figure 19f**). When used as a medial/lateral connector within the construct, the three styles of CD HORIZON® SOLERA® Connectors enable the surgeon to connect any CD HORIZON® SOLERA® rod to any CD HORIZON®


SOLERA® screw. The tightening sequence is to first tighten the connector set screw to the rod to ensure proper seating and then tighten the set screw to the pedicle screw allowing the Multi-Axial Screw head to adapt to the position of the connector post. For information on iliac fixation and connector placement, please refer to the CD HORIZON® Iliac Fixation Surgical Technique.



Figure 19f

Interbody or Vertebral Body Replacement Options

Depending on the type of surgical technique and the patient pathology, a variety of interbody or vertebral body replacement options are available. Interbody implants are to be used with autogenous bone graft and supplemental fixation for fusion in patients with degenerative disc disease. Vertebral body replacement devices are to be used with autograft or allograft and supplemental fixation to replace a collapsed, damaged, or unstable vertebral body due to tumor or trauma. For comprehensive instructions for implantation of these implants, refer to their respective surgical techniques.



*Supplemental fixation is only required for this device when the surgeon chooses to use less than three or none of the provided screws.

Rod Selection

The CD HORIZON SOLERA® Spinal System offers a spectrum of rods with different material types and lengths to facilitate construct tailoring intraoperatively (**Figure 20a**). The low profile CD HORIZON® SOLERA® screws are compatible with a variety of rod material types: Commercially Pure Titanium, Titanium Alloy, CHROMALLOY™, and CHROMALLOY™ Plus. This spectrum of rod materials and diameters allow real-time rod choices in the OR to customize the stiffness and strength of the construct to the patient needs.

Rod Options	4.75mm	5.5mm	6.0mm	5.5mm to 6.30mm by 5.8mm to 5.5mm	6.0mm to 6.30mm by 5.8mm to 6.0mm	6.30mm by 5.8mm	5.5mm to 4.75mm	6.0mm to 5.5mm
Pre-bent CHROMALLOY™	•	•	•					
Pre-bent Commercially Pure Titanium		•						
Straight Titanium Alloy	•	•	•					
Straight Commercially Pure Titanium		•	•					
Straight CHROMALLOY™	•	•						
Straight CHROMALLOY™ Plus	•	•	•					
Apex CHROMALLOY™ Plus				•	•	•		
Tapered CHROMALLOY™ Plus							•	•

Figure 20a

Pre-bent Rods

These rods range in lengths from 30mm to 120mm in 5mm increments (**Figure 20b**). The pre-bent rods reduce the steps associated with measuring, cutting, and bending straight rods during lumbar fusion surgeries. The pre-bent rods have short lines on each end for alignment during rod placement.

NOTE: Prior to implantation, break off the VERIFY® Implant Tracking Tag and retain it in the Tag Sorter so that it can be scanned at the end of the surgery.



Figure 20b

500mm and 600mm Straight Rods

There are many types of 500mm and 600mm straight rods and each rod type has different strength characteristics. The choice of the rod type depends upon the patient pathology, bone quality, and construct strength requirements as determined by the surgeon.

*The 500mm straight Titanium Alloy rod is not in the implant set and can be ordered as an extra item.

Rod Selection *continued*

Apex Rods

CHROMALOY™ Plus Apex and Tapered Rods of varying diameters are available in 500 and 600mm lengths. Apex Rods have a larger geometry in the midsection of the rod for augmented strength at the apex of a curve (**Figure 21a**). The Apex Rods have flat segments on two sides of the larger midsection and must be positioned in the screw head with the flat segments touching the sides of the screw head (**Figure 21b**). Midsection rod reduction should be considered when planning derotation and rod bending maneuvers.




Figure 21a
Transition area on both ends of Apex rod.



Figure 21b
Flat sides of Apex rod midsection seat flush with the sides of the screw head.

Rod Selection *continued*

Tapered Rods

Tapered Rods transition from a smaller diameter to a larger diameter for varied strength at the cephalad or caudal end of a construct (Figures 22a and 22b).

Tapered Rod Measurements

Length	Minor Diameter Length	Major Diameter Length
500mm	200mm at Ø4.75mm	300mm at Ø5.5mm
500mm	200mm at Ø5.5mm	300mm at Ø6.0mm
600mm	300mm at Ø4.75mm	300mm at Ø5.5mm
600mm	300mm at Ø5.5mm	300mm at Ø6.0mm



Figure 22a
Tapered rods transition from a 4.75mm to a 5.5mm or 5.5mm to a 6.0mm diameter.




Figure 22b
Diameter transition in the middle of the Tapered rod.

Rod Contouring and Placement

Once correct screw placement has been verified radiographically, measure and contour the selected rods in the sagittal and coronal planes. A rod template may be used to measure the rod length required for the construct (**Figure 23a**). A rod cutter (handheld or table top) may be used to cut the appropriate rod length.

The Titanium Alloy rods have an orientation line that serves as a reference point during contouring. Clamping the rod with Dual Action Rod Grippers at both ends helps prevent the rod from rotating during contouring (**Figure 23b**).

NOTE: Prior to implantation of the rod, break off the VERIFY® Implant Tracking Tag and retain it in the Tag Sorter so that it can be scanned at the end of the surgery.



Figure 23a



Figure 23b

Provisional Driver Assembly

The Provisional Driver assembly consists of a handle and a separate insert shaft. The instrument set contains a Provisional Driver insert shaft and a Rod Pusher insert shaft. When assembling either insert shaft, align the black laser markings on the shaft with the black line on the handle and insert the shaft (Figure 24a). When the insert shaft stops push the button on the handle and insert until it stops again (Figure 24b). Turn

the shaft 90° following the black line on the shaft (Figure 25). Push the shaft toward the handle until it locks into place. The driver is now locked in the short position (Figure 26a).

The Provisional Driver has a short and an extended position. The push button on the handle allows the driver to lock in the two different positions. The short position is used for tightening the set screw

while the extended position allows use through the Beale Rod Reducer or a Sequential Reducer.

To use the driver in the extended position, push the button and start pulling the insert shaft out from the short position. Release the button and continue pulling the shaft out until the button engages with an audible "click" (Figure 26b).




Figure 24a




Figure 24b




Figure 25

! Important

After locking the shaft into the handle, it is advised to pull back on the shaft, without depressing the button, to ensure that the shaft is securely attached to the handle.

✓ Note

This instrument requires disassembly for proper sterilization after surgery.



Figure 26a



Figure 26b

Rod Reduction

For non-hyperkyphotic deformities, place the rod on the concave side first. The contoured rod is introduced into the previously placed screws. There are several methods and instruments that facilitate fully seating the rod into the saddle of the implant.

! Important

Care should be taken with any of the following reduction methods. Improper instrument use may loosen implants or damage the residual facets and other bony anatomy.

! Important

As the Cobalt Chrome material is stiffer than Titanium material, reduction tools or manual reduction techniques must be used for rod manipulation or to seat the rod into the tulip head. Never use the set screw to reduce a spinal rod because the force applied by a set screw may not be able to fully seat the rod into the saddle of the screw.

Forceps Rocker Method

Use of the Forceps Rocker is an effective method for reducing (or seating) the rod into the implant when only a slight height difference exists between the rod and the implant saddle. To use the Forceps Rocker, grasp the sides of the implant with the rocker cam above the rod and

then lever backward over the rod (**Figure 27**). The levering action allows the rod to be fully seated into the saddle of the implant. The Dual Ended Set Screw Starter is then used to introduce the set screw (**Figure 28**).

4.75mm European Rocker* Method

The 4.75mm European Rocker* can also be used to achieve rod reduction. To use this rocker, grasp the sides of the implant with the rocker cam above the rod, squeeze the handle to secure the instrument to the implant, and then lever the instrument backward

over the rod (**Figure 29**). The levering action allows the rod to be fully seated into the implant saddle. To remove the instrument, press the blue button on the top of the instrument.




Figure 27




Figure 28




Figure 29

*Can be ordered as an extra instrument.

Rod Reduction *continued*

Beale Rod Reducer

In situations where the rod rests at the top of the implant, the Beale Rod Reducer may be used to seat the rod. The Beale Rod Reducer attaches to the four implant slots (**Figure 30**).

Once the Beale Rod Reducer is attached to the implant, squeeze the reducer handles slowly, allowing the sleeve to slide down, and seat the rod into the implant saddle.

NOTE: Prior to implantation of the set screw, break off the VERIFYI® Implant Tracking Tag and retain it in the Tag Sorter so that it can be scanned at the end of the surgery.

The set screw is then placed through the reducer tube and into the implant head with the Provisional Driver or a Dual Ended Set Screw Starter. Provisionally tighten the set screw with the Provisional Driver in the extended position (**Figure 31**).



Figure 30



Figure 31

Sequential Rod Reduction

The 4.75mm Sequential Reducer may be used to gradually seat the rod into 4.75mm implants. Insert the Inner Sleeve into the Outer Sleeve (**Figure 32**) and manually turn the wing nut clockwise until the word “LOAD” is visible in the oval window of the Sequential Reducer (**Figure 33**). Place the Sequential Reducer over the rod and rock the instrument onto the

screw head. This will allow it to “pop” into place. Then turn the wing nut until the instrument is firmly attached to the screw head. Once attached to the screw turn the wing nut clockwise until “RD” appears in the oval window and a black line is visible above the wing nut (**Figure 34**). Ensure that the rod is fully reduced using visual confirmation.



Figure 32



Figure 33



Figure 34

Sequential Rod Reduction *continued*

Attach the set screw to the Dual Ended Set Screw Starter or Provisional Driver and place it through the cannula of the reducer (**Figure 35**). Provisionally tighten the set screw and then remove the set screw starter.

To remove the Sequential Reducer from the implant turn the wing nut counterclockwise until the word "LOAD" has passed slightly below the bottom of the oval window (**Figure 36**). An audible "click" may be heard. Turning the wing nut beyond LOAD will allow the instrument to be disassembled.

✓ Note

Stop turning the handle when the "LOAD" sign passes slightly below the oval window to prevent detaching of the Inner Sleeve from the Outer Sleeve.

! Important

The Sequential Rod Reducer requires disassembly for proper sterilization after surgery.




Figure 35



Figure 36

5.5/6.0mm Lateral Translator

In situations where the rod rests medial or lateral to the top of the implant, the Lateral Translator may be used to align the implant and rod, and sequentially reduce the rod into the implant. Attach the instrument to the implant head and squeeze the handles to align the implant and rod (**Figure 37a**). Once alignment has been achieved, the spinal rod can be reduced into the implant by attaching the Translator Driver and rotating clockwise until the rod is seated (**Figure 37b**). Using the Dual Ended Set Screw Starter or Provisional Driver, a set screw can be introduced through the cannula of the Lateral Translator to provisionally tighten the set screw (**Figure 37c**). To remove the instrument from the implant, release the handle and pull the instrument up.



5.5/6.0mm SMARTLINK™ Extenders

The SMARTLINK™ Extenders may be used to gradually align and seat a 5.5mm or 6.0mm rod into 5.5/6.0mm implants. Two extenders are available, an Open and a Closed Extender. To attach either extender to a bone screw, slide the green sleeve to the unlocked position and dock the extender to the top of the screw head, aligning the implant slots to the instrument (**Figure 38a**). Once aligned, lock the Extender to the implant by pushing the green sleeve downward until the top surface of the sleeve is flush with the Extender (**Figure 38b**).

To use the Open Extender more easily with Multi-Axial Screws, engage the Extender Guide. Slide the green sleeve on the Open Extender to the unlocked position and insert the Extender Guide. The Extender Guide docks into the Multi-Axial Screw rod slot and aligns the Open Extender (**Figure 38c**). Once the Open Extender is aligned, slide the green sleeve downward into the locked position and remove the Extender Guide. Please note that the Extender Guide may not be used if the rod has already been placed.

Closed Extender




Figure 38a




Figure 38b

Open Extender with Extender Guide



Figure 38c

5.5/6.0mm SMARTLINK™ Extenders *continued*

To assemble the Reducer with either extender, align the Reducer with the internal slots of the Extender (**Figure 38**). Turn the Reducer by hand until it is engaged with the rod (**Figure 39**). Holding the top of the instrument, manipulation of the rod and implant can be done until the desired alignment is achieved. To gradually reduce the rod into the implant, the blue Driver Handle is engaged with the Reducer and manually turned clockwise (**Figure 40**).

In this seated position, remove the Driver Handle and insert a set screw to secure the rod to the bone screw. Introduction of the set screw is accomplished using the Dual Ended Set Screw Starter or the Provisional Driver. Attach the set screw to either instrument and insert it through the cannula and provisionally tighten (**Figure 41**). To remove the SMARTLINK™ Extenders from the implant, pull up on the green sleeve to disengage the tips of the instrument from the implant (**Figure 42**). Manually turning the Reducer counterclockwise will allow the Reducer to be disassembled from the Extender.



Figure 38

Figure 39

Figure 40

Figure 41

Figure 42

SMARTLINK™ Derotators

With the addition of SMARTLINK™ Derotators, the SMARTLINK™ Extenders may be used to achieve three-dimensional segmental and en bloc manipulation during complex cases. The SMARTLINK™ Derotators are designed to triangulate the SMARTLINK™ Extenders and facilitate apical vertebral body manipulation. When the SMARTLINK™ Instruments are assembled and the construct is triangulated, the Multi-Axial Screws (MAS) will mimic the axial control of a Sagittal Adjusting Screw (SAS) and a Fixed Angle Screw (FAS) allowing axial plane rotation.

To perform segmental derotation, place the Segmental Link (5485910) over bilateral SMARTLINK™ Extenders and turn the blue handle clockwise to tighten. The handle can then be used for vertebral body manipulation (**Figure 43**). If contrarotational corrective forces are desired, additional Segmental Links may be added as needed. If Multi-Axial Screws (MAS) are used on one side of the construct, place the Segmental Link over the Extender that is attached to a MAS ensuring that it passes through the round MAS hole in the Segmental Link (**Figure 44**).



Figure 43




Figure 44

SMARTLINK™ Derotators *continued*

To maintain axial control when using a Segmental Link with bilateral MAS implants, a Screw Locker* (5485912) must be used in place of the SMARTLINK™ Extender. Prior to rod placement on the convex side of the curve, place the Screw Locker onto the head of the MAS implant (Figure 45a), aligning the pusher with the rod slots, and lock in place using the Provisional Driver. The Segmental Link should then be placed over the Screw Locker and SMARTLINK™ Extender instruments ensuring that the SMARTLINK™ Extender passes through the round MAS hole (Figure 45b). Tighten the Segmental Link by turning the blue handle clockwise (Figure 45c). If the

Multi-Axial Screw does not feel completely fixed after the Segmental Link is tightened, you may reinsert the Provisional Driver and tighten further using the Segmental Link as a counter-torque. After the derotation maneuver is performed and locked in place with the concave rod, the Screw Locker may be removed and the convex rod placed.

If performing a unilateral derotation, the Interlink (5485911) should be used. Place the Interlink over multiple Extenders with the blue handle on the lateral side of the Extenders. Turn the blue handle clockwise to tighten the Interlink prior to performing derotation maneuvers (Figure 45d).



Figure 45a




Figure 45b




Figure 45c




Figure 45d

*May be ordered as an extra instrument.

SMARTLINK™ Derotators *continued*

Global derotation constructs can be assembled using a combination of Segmental Links and Interlinks or by using the Global Link (5485913), which triangulates the construct with one instrument. When using the combination of instruments, begin by placing Segmental Links over two or more levels. The Segmental Links should be pushed down the Extender shaft to accommodate the subsequent placement of the Interlink. Screw Extenders attached to MAS must pass through the MAS hole in the Segmental Link. Next, attach Interlinks on both sides of the construct to perform a Global Derotation (**Figure 46a**).

When using the Global Link, the end with the tall blue handle should be placed over the screw Extenders on the concave side first. Next, place the section with the shorter blue handle over the convex side, and tighten both handles to lock the Global Link in place. The tall blue handle may be used for vertebral manipulation. Additionally, the Lateral Handle (5485914) can be attached to either end of the Global Link, depending upon surgeon preference (**Figures 46b**). Should the sides need to move independent of each other in a maneuver, such as a rod capture, the blue handle can be loosened to allow the contralateral side to move freely.




Figure 46a




Figure 46b

Deformity Correction

The set screws are kept loose (or only locked at one end); then the concave rod is slowly straightened with the left and right Coronal Benders (**Figure 47**). Each straightening of the concave rod is performed over a pedicle screw. Several passes may be required in order for viscoelastic relaxation with subsequent curve correction to occur. Tighten the apical set screws and perform the appropriate compression or distraction. Watch the bone-to-screw interface with all correction maneuvers.

Hinged Translator

The Hinged Translator can be used in place of either a compressor or a distractor during correction maneuvers. The straight leg of the instrument will push the implant while the hinged leg engages on the rod to act as rod gripper. Pay careful attention to the bone-to-screw interface during any correction maneuver.

Prior to placing the Hinged Translator on the rod, disengage the rack so that the hinged leg and straight leg are touching each other (**Figure 48**). A left and a right translator are included in the set to facilitate the compression and distraction maneuvers around the bony anatomy. The arrow on the rack of the Hinged Translator shows the direction in which the implant will be moved.




Figure 47



Figure 48

Deformity Correction *continued*

Hinged Translator *continued*

Example for Compressing the T8-T9 Segment: Provisionally tighten the T9 set screw. Prior to squeezing the handles, place the instrument along the rod with the straight leg below and immediately against the T9 screw. (**Figure 49**). Squeeze the handles to begin compression (**Figure 50**).

Example for Distracting the T8-T9 Segment: Provisionally tighten the T8 set screw. Prior to squeezing the handles, place the instrument along the rod with the straight leg below and immediately against the T8 screw (**Figure 51**). Squeeze the handles to begin distraction (**Figure 52**).



Figure 49




Figure 50



Figure 51



Figure 52

Deformity Correction *continued*

Placing the Stabilizing Rod

Following placement of the second rod and set screws (**Figure 53**), convex compressive forces are placed on the segments using the Compressor to horizontalize the lowest instrumented vertebra and mildly compress the convexity of the deformity (**Figure 54**). It is preferred

that compression be released just prior to the set screw being broken off or with final tightening. This technique will help ensure that the implant head and rod are normalized to one another and thus allow for the rod to be fully seated in the implant head during the final

tightening step. NMEP and/or SSEP monitoring are performed to detect any potential neurologic deficits. Fixation is verified with AP and lateral x-rays to confirm spinal correction and alignment.

! Important

The Dual Ended Set Screw Starter or the Provisional Driver may be used for provisionally tightening the set screws.




Figure 53




Figure 54

Final Tightening and Decortication

Using the Counter Torque and the Self-Retaining Breakoff Driver, shear off the set screws which locks the rods into place (**Figure 55**). The break off set screw has the appropriate locking torque built into it and should not require additional tightening. Final tightening torque range is 9-10.5Nm or 80-93 in-lbs for 4.75mm implants and 10.50-12.50Nm or 92-110 in-lbs for 5.5/6.0mm implants.

If additional manipulation of the set screw is desired after the break off is achieved, the Torque Indicating Driver should be used to prevent over tightening of the set screw which could reduce the strength of the connection.

To use the 4.75mm or the 5.5/6.0mm Torque Indicating Driver, attach the Quick Connect T-Handle to the Torque

Indicating Driver and pass it through the Counter Torque and into the inner portion of the set screw (**Figure 56**). Turn the handle until the slot reaches the line on the right side of the scale to ensure the Correct Torque limit has been achieved (**Figure 57**). The posterior elements are decorticated with a burr and the bone graft is placed.

! Important

Ensure the correct Torque Indicating Driver is used for the system. A Counter Torque should be used when using the Torque Indicating Driver.



Figure 55



Figure 56



Figure 57

Bone Graft Options

Precise placement of the bone graft (autograft or allograft bone) on the decorticated surface is essential to facilitate fusion. A number of Medtronic bone graft options are available as fillers for voids or gaps that are not intrinsic to the stability of the bony structure.



MASTERGRAFT® Ceramic Scaffolds



PROGENIX® Putty



PROGENIX® Plus



MAGNIFUSE® Bone Graft
Must be placed on well decorticated surfaces.



GRAFTON® DBM Matrix Strips

CD HORIZON® X10 CROSSLINK® Plate Placement

It has been shown that a CROSSLINK® Plate provides resistance against both axial and torsional loads by converting a two rod construct into an unitized quadrilateral frame (Johnston, Ashman, Allard: Effect of Spinal Construct Stiffness on Early Fusion Mass Incorporation. *Spine* 15: 908-912, 1990). In long constructs, the CROSSLINK® Plate should be placed on the upper one-third of the construct and another one in the lower one-third of the construct.

The plates are available for 4.75mm, 5.5mm, and 6.35mm rod diameter system. To determine the appropriate CROSSLINK® Plate, use the measuring credit card or the measuring caliper (**Figure 58**).

In-line Plate Holder Method

NOTE: Prior to implantation of the CROSSLINK® Plate, break off the VERIFYI® Implant Tracking Tag and retain it in the tag sorter so that it can be scanned at the end of the surgery.

The midline nut is provisionally tightened to gain control of the CD HORIZON® X10 CROSSLINK® Multi-Span Plate. The rod set screws

are backed out such that they do not obstruct rod introduction. With the use of the In-line plate holder, the plate is gripped and positioned to capture the far rod of the two rods. The far rod set screw is provisionally tightened using the 7/32" Torque-Limiting Set Screw Driver to firmly anchor the device to the rod (**Figure 59**). Next,

loosen the midline nut to allow the multi-axial flexibility of the CD HORIZON® CROSSLINK® Plate and assemble the plate to the other rod and provisionally tighten the set screw. Retighten the midline nut to secure the overall device (**Figure 60**).



Figure 58



Figure 59



Figure 60

CD HORIZON® X10 CROSSLINK® Plate Placement *continued*

T-Bolt Implant Positioner Method

With the use of the implant positioner instruments, the appropriate CD HORIZON® X10 CROSSLINK® MULTI-SPAN Plate is selected and gripped (**Figure 61**). Ensure that both positioners fit securely onto both set screws.

The T-bolt Implant Positioners can be used to sequentially articulate the plate around the rod (**Figure 62**). If the plate cannot be precisely seated against the rod, the set screw is still too prominently extended into the claw opening. Keep the plate in the wound and abutting against the rod. By rotating the positioners, the set screw can be manipulated and slightly backed out, allowing the rod to fully seat in the claw opening. Once precise

contact has been achieved between the plate and the rod, the positioners can be used to provisionally tighten the plate to the rod. The same process is carried out for the other side of the plate. Both halves of the plate should precisely articulate with the rod before final tightening and set screw breakoff.

Remove the T-bolt Implant Positioners and provisionally tighten the midline nut using the 7/32" Torque-Limiting Set Screw Driver.



Figure 61



Figure 62

Final Tightening of the Plate

Finally tighten the midline nut to 80in-lbs (9Nm) using the Counter Torque to minimize torque transfer and the 7/32" Torque-Limiting Set Screw Driver (**Figure 63**). The 7/32" Torque-Limiting Set Screw Driver should be fully seated on the midline nut during the final tightening. The midline nut is not

a Break-off set screw. An audible "click" from the handle will confirm that the midline nut is adequately tightened to the appropriate torque.

Advance the break off set screws using the Counter Torque and the 7/32" Torque-Limiting Set Screw Driver and tighten the set screws to break off at 55-65 in-lbs

(6.2-7.3Nm) (**Figure 64**). The sheared off sections of the set screws can remain housed in the shaft of the driver until removal is convenient. To remove the sheared-off sections of the set screws from the driver, hold the handle horizontally and the broken off sections will easily fall from the oblong window in the shaft (**Figure 65**).



Figure 63



Figure 64



Figure 65

Postoperative Care and Mobilization

Prior to closure, do a final check to ensure that the set screws are symmetrically seated in the screw heads and sheared off, that the bone graft has not become dislodged during manipulation, and that a proper count of all sheared-off set screw heads is correct (**Figure 66**).

Appropriate postoperative monitoring following evaluation of the extent of the surgical procedure and the patient's overall medical status is essential. Deep vein anti-embolic treatment should be considered for all patients, along with active pulmonary toilet, fluid balance, nutritional status, and

monitoring of neurologic function. Prophylactic antibiotics may be continued for a brief duration following surgery until the wound seals. Finally, postoperative bracing may be considered for longer fusions depending upon individual surgeon preference.

A structured, progressive physical therapy program is essential to mobilize the patient in order to diminish postoperative complications and to rehabilitate the patient sufficiently for discharge. During the inpatient rehabilitation period, patients should be carefully instructed in the appropriate

methods of getting in and out of bed, stair climbing, and brace application, as well as how long to sit and various other activities of daily living. Patients who lag behind a normal recovery period proportional to the extent of their surgery should be expediently considered for transfer to a rehabilitation inpatient facility.

Finally, postoperative follow-up for a minimum of two years is crucial to assess the progression of fusion and, equally important, the patient's clinical improvement.




Figure 66



Hook Surgical Technique

Surgical Strategy

Preoperatively, any spinal surgery should be studied and a scheme of the construct defined.

Shown below are examples of some typical hook construct schemes for a T5-L2 and T2-L5 placement. These schemes, which are strictly for

illustrative purposes, are examples of how to treat various degrees of scoliosis. **Figure 67** shows a standard right thoracic curve (Lenke Type 1AN/King Type III) instrumented with hooks from T5-L2. This case can also be treated using a hybrid

construct consisting of hooks and pedicle screws. **Figure 68** shows a construct treating scoliosis from T2 to L5.




Figure 67

Hook Construct Legend	
NBH	= Narrow Blade Hook
OH	= Offset Hook
PH	= Pedicle Hook
⊗	= Pedicle Screw
WBH	= Wide Blade Hook
↗	= Up-Going Hook
↘	= Down-Going Hook
TAPH	= Total Anatomical Pedicle Hook
TATP	= Total Anatomical Transverse Process Hook
EBH	= Extended Body Hook



Figure 68

Hook Site Preparation, Options, and Insertion

The CD HORIZON® SOLERA® Spinal System is compatible with a number of top-loading hooks of different anatomic shapes and sizes (see hook implants chart, page 81). The surgeon must choose the appropriate hook based on the individual patient's anatomy, deformity degree and type, method of correction chosen, and amount of compression/distraction that will be needed to provide proper and stable purchase of the implants.

NOTE: Prior to implantation of the 4.75mm hooks, break off the VERIFY® Implant Tracking tag and retain it in the Tag Sorter so that it can be scanned at the end of the surgery.

Several different instruments can be used for hook insertion. For instance, the Straight or Lateral Implant Holder combined with the Hook Pusher (Figure 69).



Figure 69

Hook Site Preparation, Options, and Placement

Pedicle Hook

The Pedicle Hook may be used from T1 to T10. The hook blade is always cephalad (up-going) and is in the infralaminar position. The facet capsule is divided, and a portion of the inferior facet process may be removed to facilitate insertion of the hook (**Figure 70**). Once the pedicle has been clearly identified with the help of the Pedicle finder (**Figure 71**), the hook may be inserted.

If needed, a mallet can be used to impact the Hook Pusher to drive the Pedicle Hook. It is important that the Pedicle Hook is placed into the joint cavity (**Figure 72**) and is not splitting the inferior articular process (**Figure 73**).



Pedicle Hook




Figure 70



Figure 71



CORRECT
Figure 72



INCORRECT
Figure 73

Hook Site Preparation, Options, and Placement *continued*

Transverse Process Hook

This is generally a wide blade hook and is typically used in a pedicle-transverse claw construct as a caudal (down-going) as well as cephalad (up-going) hook (**Figure 74**). The Transverse Process Elevator or the wide blade Laminar Elevator may be used to separate the ligamentous attachment between the undersurface of the transverse process and the posterior arch of the rib medial to the rib-transverse joint. An Implant Holder is used to insert this hook.



Wide Blade Hook



Figure 74

Hook Site Preparation, Options, and Placement *continued*

Total Anatomical Hooks

TAH™ Total Anatomical Hooks have a small shelf designed to enhance their stability. The combination of the shelf and the close fit of the throat of these hooks demands that the angle of insertion is less vertical than required by other implants. To achieve this angle of insertion without violating the cut surface of the superior articular facet, a small amount of the adjacent inferior transverse process and lamina may need to be removed.

The TAH Pedicle Hook may be used from T1 to T10. The hook blade is always cephalad (up-going) and is in the infralaminar position. The facet capsule is divided, and a portion of the inferior facet process may be removed to facilitate insertion of the hook (**Figure 75**). Once the pedicle has been clearly identified with the help of the Pedicle Elevator (**Figure 76**), the hook may be inserted.



Total Anatomical Hook



Figure 75




Figure 76

Hook Site Preparation, Options, and Placement *continued*

TAH Transverse Process Hook

This hook is typically used in a transverse process/pedicle claw construct as a caudal (down-going) as well as cephalad (up-going) hook (**Figure 77**). The Transverse Process Elevator or the Laminar Elevator may be used to separate the ligamentous attachment between the undersurface of the transverse process and the posterior arch of the rib medial to the rib/transverse process joint. An Implant Holder is used to insert this hook.



Total Anatomical
Transverse Process Hook



Figure 77

Decortication

Once inserted, laminar hooks are not very stable prior to rod insertion. Therefore, it is recommended to remove them to decorticate.

At this point in the surgery, bilateral partial facetectomies are carried out (**Figure 78**). The intervening cartilage is denuded to allow exposure of the subchondral bone assisting in bone fusion. Decortication of the laminae,

spinous processes, and transverse processes, along with bone graft placement, will be done at the end of the surgery to avoid intraoperative bleeding. Laminar hooks are placed back into their position.



Figure 78

Rod Contouring

Once the hooks on the correction side of the deformity (concave in the thoracic area, convex in the lumbar area of the spine) are tested for fit and placement, a rod template may be used to determine the length and the curve. The correction rod is cut to the appropriate length (1cm to 2cm longer than the overall hook-to-hook length). To achieve the correct sagittal plane contour, the rod is

bent in small incremental steps using a French Bender (**Figure 79**). It is important to maintain a same plane orientation of the rod to prevent a spiral-type bend down the rod.

In the case of a reducible scoliosis, the rod is bent according to the final postoperative planned correction to obtain a nice postoperative thoracic kyphosis and lumbar lordosis.

In a case of stiff scoliosis, the rod is placed along the spine to check for proper correction, hook fit, and contouring. This type of scoliosis correction will be mainly obtained with *in situ* bending.



Figure 79

Rod Insertion

The contoured rod is placed into the top-loading implants beginning from either the upper or lower part of the construct, there is no particular rule for rod insertion. One can start with the implants in which the rod seems to best position and facilitate the continuation of the insertion (Figures 80 and 81). A Rod Holder

may be used to assist in placing the rod. Using the Dual Ended Set Screw Starter or Provisional Driver, set screws are placed into the first implants where the rod seats perfectly. The Rod Pusher may be used to push the rod down in order to place a set screw (Figure 82).

There are several methods and instruments that may be used to facilitate rod reduction and to fully seat the rod into the saddle of the implants. Refer to the Rod Reduction steps on pages 40 through 46 of the pedicle screw section of this technique for method options.




Figure 81

Figure 80



Figure 82

Deformity Correction

At this point of the surgery some of the correction has been achieved, mainly due to translation maneuvers used when inserting the rod. Further correction can be accomplished with rod rotation and/or *in situ* bending, depending on the type and stiffness of the curve, and completed with compression/distraction maneuvers.

Rod Rotation

Once the contoured rod and all of the set screws have been placed, the rod is ready to be rotated into its final position. The rotation must be done slowly in order to prevent rapid neurologic changes and/or partial pullout or hook dislodgement. The rotation is done using two Dual Action Rod Grippers (**Figure 83**). It is important to monitor the interval hooks, which

can back out during rod rotation. Several methods are proposed: use of the C-Shaped Rod Pusher, the placement of C-rings on the rod prior to rotation, placement of the Rod Gripper on the rod just below the hook to buttress it, or the use of a hook stabilizer instrument, which is available upon special ordering request.




Figure 83

Deformity Correction *continued*

Once the rotation of the rod is complete and the position of the hooks is verified, the interval hooks' set screws are provisionally tightened to prevent rod derotation. The hooks should be checked following all rotation maneuvers and the necessary adjustments made to ensure that proper placement is maintained. At this point, the rod should be fully seated into the saddle of all of the implants.

In Situ Bending

In Situ Benders may be used for correction and final adjustment of the rod in the sagittal plane. The rod is bent in small incremental steps using the two bender tips positioned near each other on the rod (**Figure 84**).




Figure 84

Compression and Distraction

Once the rod is secured in the implants, distraction and/or compression are performed to place the hooks in their final position. The Hinged Translator, Multilevel Hook Compressor, Distractor, and Provisional Driver are used to carry out these maneuvers. **(Figure 85).** Compression maneuvers are most often carried out directly on two hooks **(Figure 86).** Another option is to use the Hinged Translator for Compression. Care should be taken to ensure that the foot of either instrument is placed against the implant body and not against the set screw. It is preferred that compression be released just prior to the set screw being broken off or final tightened. This technique will help ensure that the implant head and rod are normalized to one another, and thus, allows for the rod to be fully seated in the implant head during the final tightening step. After these maneuvers are complete, the set screw is tightened with the Provisional Driver.



Figure 85




Figure 86

Stabilization and Holding Rod Placement

With the completion of the deformity correction and the seating of the correction rod, the opposite side of the construct is prepared. Measure the length for the stabilizing rod, then cut. Using the French Bender, contour the rod according to the curvature of the spine and the residual position of alignment from the correction rod. Place the contoured rod into the hooks and provisionally secure the rod with set screws (**Figure 87**). Once the rod is secured to the implants, distraction and/or compression are performed to place the hooks in their final position. Refer to the previously described instructions to ensure the appropriate steps are followed.

The spine may be decorticated to carry out the bone fusion and morselized cancellous bone placed along the decorticated spine, extending out over the transverse processes.




Figure 87

Final Tightening

When all implants are securely in place and the rod fully seated, final tightening and/or break-off of the set screws is performed.

The Counter Torque instrument is placed over the implant and the rod (**Figure 88**). The Self-Retaining Breakoff Driver is then placed through the cannulated Counter Torque. The Self-Retaining Break-Off Driver provides adequate leverage for breaking the set screw heads. The handle of the Counter Torque device should be held firmly to prevent torquing of the construct while the set screw is secured and sheared off (**Figure 89**). The broken-off part of the set screw is captured in the cannulated portion of the Self-Retaining Break-Off Driver. Following final tightening, the sheared-off portions of the set screws accumulated in the driver are removed using the T25 or T27 Obturator shaft (**Figure 90**).

! Important

The set screw should not be broken off or final tightened under compression or distraction due to possible loosening or disassembly.




Figure 88



Figure 89



Figure 90

CD HORIZON® X10 CROSSLINK® Plate Placement and Closure

Once final tightening of the set screws is completed, transverse links should be placed if possible to provide rotational stability to the construct. A framed construct resists rotational forces. Refer to the previously described instructions for placing CD HORIZON® X10 CROSSLINK® Plates in the pedicle screw section of this technique. The posterior elements should be decorticated with a burr followed by bone graft placement (**Figure 91**). Wound closure is performed in the customary manner.




Figure 91

Suggestions to Explant the CD HORIZON® SOLERA® Spinal System

Set Screws

The 4.75mm CD HORIZON® SOLERA® set screws (plugs) may be removed using the T25 Obturator and the Self-retaining Breakoff Driver. The T25 Obturator is inserted into the working end of the Self-retaining Break-off Driver, so that the knurled portion of the T25 Obturator is flush with the driver. Insert the obturator tip through the Counter Torque, which should be seated on the screw and into the plug, turning counterclockwise until the plug has been removed.

When removing the 5.5/6.0mm CD HORIZON® SOLERA® set screws (plug), use the T27 Obturator and following the same steps as listed above.

Pedicle Screws

The pedicle screws may be removed using either the Ball-ended T25 Bone Screw Removal Driver or the Self-retaining Screwdriver in conjunction with the Quick Connect Handle. First, attach the Quick Connect Handle to the modular end of the driver. Next, fully engage the T25 end of the driver into the screw head; then, if utilizing the Multi-Axial Screwdriver, thread the instrument sleeve into the screw head. Turn counterclockwise until the pedicle screws have been removed.

Hooks

The hooks may be removed using the Self-retaining Implant Holder. Attach the Self-retaining Implant Holder to the implant and remove the hook.

CD HORIZON® CROSSLINK® Plates

If removal of a CD HORIZON® X10 CROSSLINK® MULTI-SPAN Plate is necessary, place the 7/32" Torque-Limiting Set Screw Driver over the midline nut and turn counterclockwise to loosen. Place the 3.0mm Hex Removal Driver into a standard Medtronic Quick Connect Handle. Place the tip of the 3.0mm internal hex driver into the set screw and confirm that the 3.0mm tip is completely inserted and seated in the set screw so that the tip does not strip the hex. Turn the driver counterclockwise to loosen the set screw from the rod.

Cobalt Chrome/Titanium 5.5/6.0mm Multi-Axial Screws

May be ordered as extras

Catalog Number	Description	
55840006500	6.5mm x 100mm	
55840007500	7.5mm x 100mm	
55840008500	8.5mm x 100mm	
55840009500	9.5mm x 100mm	
55840001540	10.5mm x 40mm	
55840001545	10.5mm x 45mm	
55840001550	10.5mm x 50mm	
55840001555	10.5mm x 55mm	
55840001560	10.5mm x 60mm	
55840001565	10.5mm x 65mm	
55840001570	10.5mm x 70mm	
55840001575	10.5mm x 75mm	
55840001580	10.5mm x 80mm	
55840001585	10.5mm x 85mm	
55840001590	10.5mm x 90mm	

4.0mm	4.5mm	5.0mm	5.5mm	6.0mm	6.5mm	7.5mm	8.5mm	9.5mm	10.5mm

Product Ordering Information *continued*

5.5/6.0mm Sagittal Adjusting Screw



Catalog Number	Description	
55811004025	4.0mm X 25mm	○
55811004030	4.0mm X 30mm	●
55811004035	4.0mm X 35mm	○
55811004040	4.0mm X 40mm	●
55811004045	4.0mm X 45mm	○
55811004520	4.5mm x 20mm	●
55811004525	4.5mm x 25mm	●
55811004530	4.5mm x 30mm	●
55811004535	4.5mm x 35mm	●
55811004540	4.5mm x 40mm	●
55811004545	4.5mm x 45mm	●
55811005020	5.0mm x 20mm	●
55811005025	5.0mm x 25mm	●
55811005030	5.0mm x 30mm	●
55811005035	5.0mm x 35mm	●
55811005040	5.0mm x 40mm	●
55811005045	5.0mm x 45mm	●
55811005530	5.5mm x 30mm	●
55811005535	5.5mm x 35mm	●
55811005540	5.5mm x 40mm	●
55811005545	5.5mm x 45mm	●
55811005550	5.5mm x 50mm	●
55811006530	6.5mm x 30mm	●
55811006535	6.5mm x 35mm	●
55811006540	6.5mm x 40mm	●
55811006545	6.5mm x 45mm	●
55811006550	6.5mm x 50mm	●
55811006555	6.5mm x 55mm	●
55811007535	7.5mm x 35mm	●
55811007540	7.5mm x 40mm	●
55811007545	7.5mm x 45mm	●
55811007550	7.5mm x 50mm	●
55811007555	7.5mm x 55mm	●

Titanium Set Screws



Catalog Number	Description
5440030	4.75mm Break-off Set Screws
5540030	5.5/6.0mm Break-off Set Screws
5440130	4.75mm Non-Break-off Set Screws
5540130	5.5/6.0mm Non-Break-off Set Screws
5440230	4.75mm Reduction Set Screws
5540230	5.5/6.0mm Reduction Set Screws

Titanium Fixed Angle Screws



4.75mm	5.5mm/6.0mm	Description	
	55410004520	4.5mm x 20mm	●
54410004525	55410004525	4.5mm x 25mm	●
54410004530	55410004530	4.5mm x 30mm	●
54410004535	55410004535	4.5mm x 35mm	●
54410004540	55410004540	4.5mm x 40mm	●
54410004545	55410004545	4.5mm x 45mm	●
	55410005020	5.0mm x 20mm	●
54410005025	55410005025	5.0mm x 25mm	●
54410005030	55410005030	5.0mm x 30mm	●
54410005035	55410005035	5.0mm x 35mm	●
54410005040	55410005040	5.0mm x 40mm	●
54410005045	55410005045	5.0mm x 45mm	●
54410005530	55410005530	5.5mm x 30mm	●
54410005535	55410005535	5.5mm x 35mm	●
54410005540	55410005540	5.5mm x 40mm	●
54410005545	55410005545	5.5mm x 45mm	●
54410005550	55410005550	5.5mm x 50mm	●
54410006530	55410006530	6.5mm x 30mm	●
54410006535	55410006535	6.5mm x 35mm	●
54410006540	55410006540	6.5mm x 40mm	●
54410006545	55410006545	6.5mm x 45mm	●
54410006550	55410006550	6.5mm x 50mm	●
54410006555	55410006555	6.5mm x 55mm	●
54410007535	55410007535	7.5mm x 35mm	●
54410007540	55410007540	7.5mm x 40mm	●
54410007545	55410007545	7.5mm x 45mm	●
54410007550	55410007550	7.5mm x 50mm	●
54410007555	55410007555	7.5mm x 55mm	●

Titanium 5.5/6.0mm Fixed Angle Screws

May be ordered as extras

Catalog Number	Description	
55410004020	4.0mm x 20mm	○
55410004025	4.0mm x 25mm	●
55410004030	4.0mm x 30mm	○
55410004035	4.0mm x 35mm	●
55410004040	4.0mm x 40mm	○
55410004045	4.0mm x 45mm	●
55410005520	5.5mm x 20mm	●
55410005525	5.5mm x 25mm	●
55410006025	6.0mm x 25mm	○
55410006030	6.0mm x 30mm	○
55410006035	6.0mm x 35mm	○
55410006040	6.0mm x 40mm	○
55410006045	6.0mm x 45mm	○
55410006050	6.0mm x 50mm	○
55410006520	6.5mm x 20mm	●
55410006525	6.5mm x 25mm	●

Product Ordering Information *continued*

Multi-Axial Reduction Screws (MARS)

4.75mm	5.5mm/6.0mm	Description	
54890004525	55890004525	4.5mm x 25mm	●
54890004530	55890004530	4.5mm x 30mm	●
54890004535	55890004535	4.5mm x 35mm	●
54890004540	55890004540	4.5mm x 40mm	●
54890005525	55890005525	5.5mm x 25mm	●
54890005530	55890005530	5.5mm x 30mm	●
54890005535	55890005535	5.5mm x 35mm	●
54890005540	55890005540	5.5mm x 40mm	●
54890005545	55890005545	5.5mm x 45mm	●
54890005550	55890005550	5.5mm x 50mm	●
54890006530	55890006530	6.5mm x 30mm	●
54890006535	55890006535	6.5mm x 35mm	●
54890006540	55890006540	6.5mm x 40mm	●
54890006545	55890006545	6.5mm x 45mm	●
54890006550	55890006550	6.5mm x 50mm	●
54890006555	55890006555	6.5mm x 55mm	●
54890007535	55890007535	7.5mm x 35mm	●
54890007540	55890007540	7.5mm x 40mm	●
54890007545	55890007545	7.5mm x 45mm	●
54890007550	55890007550	7.5mm x 50mm	●
54890007555	55890007555	7.5mm x 55mm	●

Hydroxyapatite (HA) Coated Multi-Axial Screws

4.75mm	5.5mm/6.0mm	Description	
54740105535	55740105535	5.5mm x 35mm	●
54740105540	55740105540	5.5mm x 40mm	●
54740105545	55740105545	5.5mm x 45mm	●
54740105550	55740105550	5.5mm x 50mm	●
54740106540	55740106540	6.5mm x 40mm	●
54740106545	55740106545	6.5mm x 45mm	●
54740106550	55740106550	6.5mm x 50mm	●
54740106555	55740106555	6.5mm x 55mm	●
54740107540	55740107540	7.5mm x 40mm	●
54740107545	55740107545	7.5mm x 45mm	●
54740107550	55740107550	7.5mm x 50mm	●
54740107555	55740107555	7.5mm x 55mm	●
54740108535	55740108535	8.5mm x 35mm	●
54740108540	55740108540	8.5mm x 40mm	●
54740108545	55740108545	8.5mm x 45mm	●
54740108550	55740108550	8.5mm x 50mm	●

4.75mm Cannulated Screws

Catalog Number	Description	
54840015535	5.5mm x 35mm	●
54840015540	5.5mm x 40mm	●
54840015545	5.5mm x 45mm	●
54840015550	5.5mm x 50mm	●
54840016535	6.5mm x 35mm	●
54840016540	6.5mm x 40mm	●
54840016545	6.5mm x 45mm	●
54840016550	6.5mm x 50mm	●
54840016555	6.5mm x 55mm	●
54840017530	7.5mm x 30mm	●
54840017535	7.5mm x 35mm	●
54840017540	7.5mm x 40mm	●
54840017545	7.5mm x 45mm	●

5.5/6.0mm Cannulated Sagittal Adjusting Screw

Catalog Number	Description	
55811014530	4.5mm x 30mm	●
55811014535	4.5mm x 35mm	●
55811014540	4.5mm x 40mm	●
55811014545	4.5mm x 45mm	●
55811015030	5.0mm x 30mm	●
55811015035	5.0mm x 35mm	●
55811015040	5.0mm x 40mm	●
55811015045	5.0mm x 45mm	●
55811015530	5.5mm x 30mm	●
55811015535	5.5mm x 35mm	●
55811015540	5.5mm x 40mm	●
55811015545	5.5mm x 45mm	●
55811015550	5.5mm x 50mm	●
55811016530	6.5mm x 30mm	●
55811016535	6.5mm x 35mm	●
55811016540	6.5mm x 40mm	●
55811016545	6.5mm x 45mm	●
55811016550	6.5mm x 50mm	●
55811016555	6.5mm x 55mm	●
55811017535	7.5mm x 35mm	●
55811017540	7.5mm x 40mm	●
55811017545	7.5mm x 45mm	●
55811017550	7.5mm x 50mm	●
55811017555	7.5mm x 55mm	●

Product Ordering Information *continued*

5.5/6.0mm Multi-Axial Reduction Screws (MARS)

May be ordered as extras

Catalog Number	Description	
55890004545	4.5mm x 45mm	●
55890004550	4.5mm x 50mm	●
55890004555	4.5mm x 55mm	●
55890005555	5.5mm x 55mm	●
55890005560	5.5mm x 60mm	●
55890006525	6.5mm x 25mm	●
55890006560	6.5mm x 60mm	●
55890006565	6.5mm x 65mm	●
55890007555	7.5mm x 55mm	●
55890007560	7.5mm x 60mm	●
55890007565	7.5mm x 65mm	●
55890008535	8.5mm x 35mm	●
55890008540	8.5mm x 40mm	●
55890008545	8.5mm x 45mm	●
55890008550	8.5mm x 50mm	●
55890008555	8.5mm x 55mm	●

4.75mm Multi-Axial Awl-Tap Screw (ATS™)



Catalog Number	Description	
54840044525	4.5mm x 25mm	●
54840044530	4.5mm x 30mm	●
54840044535	4.5mm x 35mm	●
54840044540	4.5mm x 40mm	●
54840045530	5.5mm x 30mm	●
54840045535	5.5mm x 35mm	●
54840045540	5.5mm x 40mm	●
54840045545	5.5mm x 45mm	●
54840045550	5.5mm x 50mm	●
54840046530	6.5mm x 30mm	●
54840046535	6.5mm x 35mm	●
54840046540	6.5mm x 40mm	●
54840046545	6.5mm x 45mm	●
54840046550	6.5mm x 50mm	●
54840047535	7.5mm x 35mm	●
54840047540	7.5mm x 40mm	●
54840047545	7.5mm x 45mm	●
54840047550	7.5mm x 50mm	●
54840047555	7.5mm x 55mm	●

CD HORIZON® SOLERA® SEXTANT® Set Screw



Catalog Number	Description
5440430	4.75mm Internal Hex Set Screws

CD HORIZON® SOLERA® SEXTANT® 4.75mm Rods



Cobalt Chrome Pre-bent Rod

Catalog Number	Description
1475005030	30mm Pre-bent Cobalt Chrome
1475005035	35mm Pre-bent Cobalt Chrome
1475005040	40mm Pre-bent Cobalt Chrome
1475005045	45mm Pre-bent Cobalt Chrome
1475005050	50mm Pre-bent Cobalt Chrome
1475005055	55mm Pre-bent Cobalt Chrome
1475005060	60mm Pre-bent Cobalt Chrome
1475005065	65mm Pre-bent Cobalt Chrome
1475005070	70mm Pre-bent Cobalt Chrome
1475005075	75mm Pre-bent Cobalt Chrome
1475005080	80mm Pre-bent Cobalt Chrome
1475005085	85mm Pre-bent Cobalt Chrome
1475005090	90mm Pre-bent Cobalt Chrome

Product Ordering Information *continued*

CHROMALLOY® Plus Apex Rods



CHROMALLOY™ Plus Apex Rod

Catalog Number	Description
1556009500*	5.5mm to 6.30mm by 5.8mm to 5.5mm CHROMALLOY™ Plus, 500mm Length
1556009600*	5.5mm to 6.30mm by 5.8mm to 5.5mm CHROMALLOY™ Plus, 600mm Length
1606009500*	6.0mm to 6.30mm by 5.8mm to 6.0mm CHROMALLOY™ Plus, 500mm Length
1606009600*	6.0mm to 6.30mm by 5.8mm to 6.0mm CHROMALLOY™ Plus, 600mm Length
1636009600*	6.3mm by 5.80mm CHROMALLOY™ Plus, 600mm Length

CHROMALLOY® Plus Tapered Rods



CHROMALLOY™ Plus Tapered Rod

Catalog Number	Description
1546000500*	5.5mm to 4.75mm CHROMALLOY™ Plus, 500mm Length
1546000600*	5.5mm to 4.75mm CHROMALLOY™ Plus, 600mm Length
1656000500*	6.0mm to 5.5mm CHROMALLOY™ Plus, 500mm Length
1656000600*	6.0mm to 5.5mm CHROMALLOY™ Plus, 600mm Length

4.75mm/5.5mm/6.0mm Rods



Commercially Pure Titanium Pre-bent Rod



CHROMALLOY™ Pre-bent Rod



Commercially Pure Titanium Rod



Titanium Alloy Rod



CHROMALLOY™ Rod



CHROMALLOY™ Plus Rod

4.75mm	5.5mm	6.0mm	Description
	1553201030		30mm Pre-bent Commercially Pure Titanium
	1553201035		35mm Pre-bent Commercially Pure Titanium
	1553201040		40mm Pre-bent Commercially Pure Titanium
	1553201045		45mm Pre-bent Commercially Pure Titanium
	1553201050		50mm Pre-bent Commercially Pure Titanium
	1553201055		55mm Pre-bent Commercially Pure Titanium
	1553201060		60mm Pre-bent Commercially Pure Titanium
	1553201070		70mm Pre-bent Commercially Pure Titanium
	1553201080		80mm Pre-bent Commercially Pure Titanium
	1553201090		90mm Pre-bent Commercially Pure Titanium
	1553201100		100mm Pre-bent Commercially Pure Titanium
	1553201110		110mm Pre-bent Commercially Pure Titanium
	1553201120		120mm Pre-bent Commercially Pure Titanium
1475501030	1555501030	1605501030	30mm Pre-bent CHROMALLOY™
1475501035	1555501035	1605501035	35mm Pre-bent CHROMALLOY™
1475501040	1555501040		40mm Pre-bent CHROMALLOY™
1475501045	1555501045		45mm Pre-bent CHROMALLOY™
1475501050	1555501050		50mm Pre-bent CHROMALLOY™
1475501055	1555501055		55mm Pre-bent CHROMALLOY™
1475501060	1555501060		60mm Pre-bent CHROMALLOY™
1475501070	1555501070		70mm Pre-bent CHROMALLOY™
1475501080	1555501080		80mm Pre-bent CHROMALLOY™
1475501090	1555501090		90mm Pre-bent CHROMALLOY™
1475501100	1555501100		100mm Pre-bent CHROMALLOY™
1475501110	1555501110		110mm Pre-bent CHROMALLOY™
1475501120	1555501120	1605501120	120 mm Pre-bent CHROMALLOY™
	1553200500	1603200500*	500mm Straight Commercially Pure Titanium
1474000500	1554200500	1604200500*	500mm Straight Titanium Alloy
1475000500	1555200500	1605000500*	500mm Straight CHROMALLOY™
1476000500	1556200500	1606200500	500mm Straight CHROMALLOY™ Plus
1476200600		1606000600	600mm Straight CHROMALLOY™ Plus
		1606000700	700mm Straight CHROMALLOY™ Plus

*May be ordered as extras

Product Ordering Information *continued*

Titanium Hooks

	4.75mm	5.5/6.0mm	Description	
	5441101	5541101*	Pedicule Hook, Extra Small	
	5441102	5541102	Pedicule Hook, Small	
	5441103	5541103	Pedicule Hook, Medium	
	5441104	5541104	Pedicule Hook, Large	
	5441112	5541112	Wide Blade Hook, Small	
	5441113	5541113	Wide Blade Hook, Medium	
	5441114	5541114	Wide Blade Hook, Large	
	5441122	5541122	Narrow Blade Hook, Small	
	5441123	5541123	Narrow Blade Hook, Medium	
	5441124	5541124	Narrow Blade Hook, Large	
	5441133	5541133	Ramped, Wide Blade, Medium	
	5441142	5541142	Ramped, Narrow Blade, Small	
	5441143	5541143	Ramped, Narrow Blade, Medium	
	5441172	5541172	Extended Body Hook, Small	
	5441173	5541173	Extended Body Hook, Medium	
	5441174	5541174	Extended Body Hook, Large	
		5541188*	Right Thoracic Angled Blade, Small	
		5541189*	Left Thoracic Angled Blade, Small	
	5441196		Right Offset Hook, Small	
	5441197		Left Offset Hook, Small	
	5441198	5541198	Right Offset Hook, Large	
	5441199	5541199	Left Offset Hook, Large	
	5441302	5541302	Total Anatomical Pedicle Hook, Small	
	5441303	5541303	Total Anatomical Pedicle Hook, Medium	
	5441342	5541342	Total Anatomical Transverse Process Hook, Small	
	5441343	5541343	Total Anatomical Transverse Process Hook, Medium	
		5541190*	Right Thoracic Angled Blade, Medium	
		5541191*	Left Thoracic Angled Blade, Medium	
		5541312*	Total Anatomical Supralaminar, Small	
		5541313*	Total Anatomical Supralaminar, Medium	
		5541322*	Total Anatomical Infralaminar, Small	
		5541323*	Total Anatomical Infralaminar, Medium	

Hook Color-coding Size Reference

Extra Small	Small	Medium	Large
			

*May be ordered as extras

Product Ordering Information *continued*

CD HORIZON® X10 CROSSLINK® Plates, Titanium

	4.75mm	5.5mm	6.35mm	Description
	5442016	8115516	8116416	16mm Fixed
	5442019	8115519	8116419	19mm Fixed
	5442022	8115522	8116422	22mm Fixed
	5442025	8115525	8116425	25mm Fixed
	5442028			28mm Fixed
	5442031			31mm Fixed
	5442034			34mm Fixed
	5442037			37mm Fixed
		8115528	8116428	28mm-30mm MULTISPAN
	5442130			30mm-32mm MULTISPAN
		8115530	8116430	30mm-34mm MULTISPAN
	5442132			32mm-37mm MULTISPAN
		8115534	8116434	34mm-36 mm MULTISPAN
	5442136			36mm-38mm MULTISPAN
		8115536	8116436	36mm-39mm MULTISPAN
	5442138			38mm-42mm MULTISPAN
		8115539	8116439	39mm-45mm MULTISPAN
	5442141			41mm-48mm MULTISPAN
	5442144			44mm-54mm MULTISPAN
		8115545	8116445	45mm-58mm MULTISPAN
	5442147			47mm-60mm MULTISPAN
	5442152			52mm-70mm MULTISPAN
		8115558	8116458	58mm-80mm MULTISPAN
	5442160			60mm-86mm MULTISPAN
	8110855	8110855	8110855	Set Screw

Product Ordering Information *continued*

Titanium Closed Multi-Axial Iliac Screws*



Catalog Number	Description	Catalog Number	Description
70465540	5.5mm × 40mm	70467560	7.5mm × 60mm
70465550	5.5mm × 50mm	70467570	7.5mm × 70mm
70465560	5.5mm × 60mm	70467580	7.5mm × 80mm
70466550	6.5mm × 50mm	70468570	8.5mm × 70mm
70466560	6.5mm × 60mm	70468580	8.5mm × 80mm
70466570	6.5mm × 70mm	70468590	8.5mm × 90mm


Titanium Iliac Set Screws




Catalog Number	Description
7049855	Hex Break-off Set Screws

Product Ordering Information *continued*


Side-Loading Lateral Connectors




Catalog Number	Description
4.75mm to 4.75mm	
5444220	4.75mm Side-Loading Lateral Connector, 4.75mm Post, 20mm Length
5444225	4.75mm Side-Loading Lateral Connector, 4.75mm Post, 25mm Length
5444235	4.75mm Side-Loading Lateral Connector, 4.75mm Post, 35mm Length
5444270*	4.75mm Side-Loading Lateral Connector, 4.75mm Post, 70mm Length



Catalog Number	Description
4.75mm to 5.5mm	
5454220	4.75mm Side-Loading Lateral Connector, 5.5mm Post, 20mm Length
5454225	4.75mm Side-Loading Lateral Connector, 5.5mm Post, 25mm Length
5454235	4.75mm Side-Loading Lateral Connector, 5.5mm Post, 35mm Length
5454270*	4.75mm Side-Loading Lateral Connector, 5.5mm Post, 70mm Length



Catalog Number	Description
5.5/6.0mm to 5.5mm	
5554220	5.5/6.0mm Side-Loading Lateral Connector, 5.5mm Post, 20mm Length
5554225	5.5/6.0mm Side-Loading Lateral Connector, 5.5mm Post, 25mm Length
5554235	5.5/6.0mm Side-Loading Lateral Connector, 5.5mm Post, 35mm Length
5554270*	5.5/6.0mm Side-Loading Lateral Connector, 5.5mm Post, 70mm Length



Catalog Number	Description
5.5/6.0mm to 6.0mm	
5564220	5.5/6.0mm Side-Loading Lateral Connector, 6.0mm Post, 20mm Length
5564225	5.5/6.0mm Side-Loading Lateral Connector, 6.0mm Post, 25mm Length
5564235	5.5/6.0mm Side-Loading Lateral Connector, 6.0mm Post, 35mm Length
5564270*	5.5/6.0mm Side-Loading Lateral Connector, 6.0mm Post, 70mm Length


Color-coding Size Reference

4.75mm	5.5/6.0mm
	


*May be ordered as extras.

Product Ordering Information *continued*


Closed Lateral Connectors



Catalog Number	Description
4.75mm to 4.75mm	
5444120	4.75mm Closed Lateral Connector, 4.75mm Post, 20mm Length
5444125	4.75mm Closed Lateral Connector, 4.75mm Post, 25mm Length
5444135	4.75mm Closed Lateral Connector, 4.75mm Post, 35mm Length
5444170*	4.75mm Closed Lateral Connector, 4.75mm Post, 70mm Length




4.75mm to 5.5mm	
5454120	4.75mm Closed Lateral Connector, 5.5mm Post, 20mm Length
5454125	4.75mm Closed Lateral Connector, 5.5mm Post, 25mm Length
5454135	4.75mm Closed Lateral Connector, 5.5mm Post, 35mm Length
5454170*	4.75mm Closed Lateral Connector, 5.5mm Post, 70mm Length




5.5/6.0mm to 6.0mm	
5564120	5.5/6.0mm Closed Lateral Connector, 6.0mm Post, 20mm Length
5564125	5.5/6.0mm Closed Lateral Connector, 6.0mm Post, 25mm Length
5564135	5.5/6.0mm Closed Lateral Connector, 6.0mm Post, 35mm Length
5564170*	5.5/6.0mm Closed Lateral Connector, 6.0mm Post, 70mm Length


Top-Loading Lateral Connectors



Catalog Number	Description
4.75mm to 4.75mm	
5444320	4.75mm Top-Loading Lateral Connector, 4.75mm Post, 20mm Length
5444325	4.75mm Top-Loading Lateral Connector, 4.75mm Post, 25mm Length
5444335	4.75mm Top-Loading Lateral Connector, 4.75mm Post, 35mm Length
5444370*	4.75mm Top-Loading Lateral Connector, 4.75mm Post, 70mm Length



4.75mm to 6.0mm	
5464320	4.75mm Top-Loading Lateral Connector, 6.0mm Post, 20mm Length
5464325	4.75mm Top-Loading Lateral Connector, 6.0mm Post, 25mm Length
5464335	4.75mm Top-Loading Lateral Connector, 6.0mm Post, 35mm Length
5464370*	4.75mm Top-Loading Lateral Connector, 6.0mm Post, 70mm Length



5.5/6.0mm to 6.0mm	
5564320	5.5/6.0mm Top-Loading Lateral Connector, 6.0mm Post, 20mm Length
5564325	5.5/6.0mm Top-Loading Lateral Connector, 6.0mm Post, 25mm Length
5564335	5.5/6.0mm Top-Loading Lateral Connector, 6.0mm Post, 35mm Length
5564370*	5.5/6.0mm Top-Loading Lateral Connector, 6.0mm Post, 70mm Length

Lateral Connector Set Screws for Closed and Side-Loading Connectors



Catalog Number	Description
779170005	Lateral Connector Set Screws

Lateral Connector Set Screws for Top-Loading Connectors



Catalog Number	Description
5440030	Lateral Connector Set Screws, 4.75mm
5540030	Lateral Connector Set Screws, 5.5/6.0mm

*May be ordered as extras.

Important Information on the CD HORIZON® Spinal System

PURPOSE

The CD HORIZON® Spinal System is intended to help provide immobilization and stabilization of spinal segments as an adjunct to fusion of the thoracic, lumbar, and/or sacral spine.

DESCRIPTION

The CD HORIZON® Spinal System consists of a variety of shapes and sizes of rods, hooks, screws, CROSSLINK® Plates, staples and connecting components, as well as implant components from other Medtronic spinal systems, which can be rigidly locked into a variety of configurations, with each construct being tailor-made for the individual case.

A subset of CD HORIZON® Spinal System components may be used for posterior pedicle screw fixation in pediatric cases. These constructs may be comprised of a variety of shapes and sizes of rods (ranging in diameter from 3.5mm to 6.35mm), hooks, screws, CROSSLINK® Plates, and connecting components. Similarly to the CD HORIZON® implants used in adult cases, these components can be rigidly locked into a variety of configurations, with each construct being tailor-made for the individual case.

Certain components within the CD HORIZON® Spinal System are specifically excluded for use in pediatric patients. These include PEEK rods, Shape Memory Alloy Staples, SPIRE™ Plates and DYNALOK bolts. All screws used in pediatric cases are only cleared for use via a posterior approach. All of the components used in pediatric cases are fabricated from medical grade stainless steel, medical grade titanium, titanium alloy, and medical grade cobalt-chromium-molybdenum alloy.

Certain implant components from other Medtronic spinal systems can be used with the CD HORIZON® Spinal System in non-pediatric cases. These components include TSRH® rods, hooks, screws, plates, CROSSLINK® plates, connectors, staples, washers, GDLH rods, hooks, connectors and CROSSLINK® bar and connectors; LIBERTY rods and screws; DYNALOK PLUS and DYNALOK CLASSIC bolts along with rod/bolt connectors; and Medtronic Multi-Axial rods and screws. Please note that certain components are specifically designed to connect to specific rod diameters, while other components can connect to multiple rod diameters. Care should be taken so that the correct components are used in the spinal construct.

CD HORIZON® hooks are intended for posterior use only. CD HORIZON® staples and CD HORIZON® ECLIPSE™ rods and associated screws are intended for anterior use only. However, for patients of smaller stature and pediatric patients, CD HORIZON® 4.5mm rods and associated components may be used posteriorly.

The CD HORIZON® Spinal System implant components are fabricated from medical grade stainless steel, medical grade titanium, titanium alloy, medical grade cobalt-chromium-molybdenum alloy, or medical grade PEEK OPTIMA-LT1. Certain CD HORIZON® Spinal System components may be coated with hydroxyapatite. No warranties expressed or implied are made. Implied warranties of merchantability and fitness for a particular purpose or use are specifically excluded. See the MDT Catalog for further information about warranties and limitations of liability.

Never use stainless steel and titanium implant components in the same construct.

Medical grade titanium, titanium alloy, and/or medical grade cobalt-chromium-molybdenum alloy may be used together. Never use titanium, titanium alloy, and/or medical grade cobalt-chromium-molybdenum alloy with stainless steel in the same construct.

The CD HORIZON® Spinal System also includes anterior staples made of Shape Memory Alloy (Nitinol – NiTi). Shape Memory Alloy is compatible with titanium, titanium alloy, and cobalt-chromium-molybdenum alloy. Do not use with stainless steel. **These staples are not to be used in pediatric patients.**

PEEK OPTIMA-LT1 implants may be used with stainless steel, titanium, or cobalt-chromium-molybdenum alloy implants. **CD HORIZON® PEEK Rods are not to be used with CROSSLINK® Plates or in pediatric patients.**

To achieve best results, do not use any of the CD HORIZON® Spinal System implant components with components from any other system or manufacturer unless specifically allowed to do so in this or another Medtronic document. As with all orthopaedic and neurosurgical implants, none of the CD HORIZON® Spinal System components should ever be reused under any circumstances.

INDICATIONS

The CD HORIZON® Spinal System with or without SEXTANT® instrumentation is intended for posterior, non-cervical fixation as an adjunct to fusion for the following indications: degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies), spondylolisthesis, trauma (i.e., fracture or dislocation), spinal stenosis, curvatures (i.e., scoliosis, kyphosis, or lordosis), tumor, pseudarthrosis, and/or failed previous fusion.

Except for hooks, when used as an anterolateral thoracic/lumbar system, the CD HORIZON® Spinal System may also be used for the same indications as an adjunct to fusion.

With the exception of degenerative disc disease, the CD HORIZON® LEGACY™ 3.5mm rods and the CD HORIZON® Spinal System PEEK rods and associated components may be used for the aforementioned indications in skeletally mature patients as an adjunct to fusion. The 3.5mm rods may be used for the specific pediatric indications noted below.

When used for posterior non-cervical pedicle screw fixation in pediatric patients, the CD HORIZON® Spinal System implants are indicated as an adjunct to fusion to treat progressive spinal deformities (i.e., scoliosis, kyphosis, or lordosis) including idiopathic scoliosis, neuromuscular scoliosis, and congenital scoliosis. Additionally, the CD HORIZON® Spinal System is intended to treat pediatric patients diagnosed with the following conditions: spondylolisthesis/spondylolysis, fracture caused by tumor and/or trauma, pseudarthrosis, and/or failed previous fusion. These devices are to be used with autograft and/or allograft. Pediatric pedicle screw fixation is limited to a posterior approach.

The CD HORIZON® SPIRE™ Plate is a posterior, single-level, non-pedicle supplemental fixation device intended for use in the non-cervical spine (T1-S1) as an adjunct to fusion in skeletally mature patients. It is intended for plate fixation/attachment to spinous processes for the purpose of achieving supplemental fixation in the following conditions: degenerative disc disease (as previously defined), spondylolisthesis, trauma, and/or tumor.

In order to achieve additional levels of fixation, the CD HORIZON® Spinal System rods may be connected to the VERTEX® Reconstruction System with the VERTEX® rod connector. Refer to the VERTEX® Reconstruction System Package Insert for a list of the VERTEX® indications of use.

CONTRAINDICATIONS

Contraindications include, but are not limited to:

- Active infectious process or significant risk of infection (immunocompromise).
- Signs of local inflammation.
- Fever or leukocytosis.
- Morbid obesity.
- Pregnancy.
- Mental illness.
- Grossly distorted anatomy caused by congenital abnormalities.
- Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of congenital abnormalities, elevation of sedimentation rate unexplained by other diseases, elevation of white blood count (WBC), or a marked left shift in the WBC differential count.
- Suspected or documented metal allergy or intolerance.
- Any case not needing a bone graft and fusion.
- Any case where the implant components selected for use would be too large or too small to achieve a successful result.
- Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality.
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance.
- The CD HORIZON® SPIRE™ Plate and the CD HORIZON® PEEK Rods are specifically contraindicated for use in pediatric patients.
- Any patient unwilling to follow postoperative instructions.
- Any case not described in the indications.

NOTA BENE: Although not absolute contraindications, conditions to be considered as potential factors for not using this device include:

• **Severe bone resorption.**

• **Osteomalacia**

• **Severe osteoporosis.**

POTENTIAL ADVERSE EVENTS

All of the possible adverse events associated with spinal fusion surgery without instrumentation are possible. With instrumentation, a listing of potential adverse events includes, but is not limited to:

- Early or late loosening of any or all of the components.
- Disassembly, bending, and/or breakage of any or all of the components.
- Foreign body (allergic) reaction to implants, debris, corrosion products (from crevice, fretting, and/or general corrosion), including metallosis, staining, tumor formation, and/or autoimmune disease.
- Pressure on the skin from component parts in patients with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, fibrosis, neurosis, and/or pain.
- Bursitis.
- Tissue or nerve damage caused by improper positioning and placement of implants or instruments.
- Post-operative change in spinal curvature, loss of correction, height, and/or reduction.
- Infection.
- Dural tears, pseudomeningocele, fistula, persistent CSF leakage, meningitis.
- Loss of neurological function (e.g., sensory and/or motor), including paralysis (complete or incomplete), dysesthesias, hyperesthesia, anesthesia, paresthesia, appearance of radiculopathy, and/or the development or continuation of pain, numbness, neuroma, spasms, sensory loss, tingling sensation, and/or visual deficits.
- Cauda equina syndrome, neuropathy, neurological deficits (transient or permanent), paraplegia, paraparesis, reflex deficits, irritation, arachnoiditis, and/or muscle loss.
- Urinary retention or loss of bladder control or other types of urological system compromise.
- Scar formation possibly causing neurological compromise or compression around nerves and/or pain.
- Fracture, microfracture, resorption, damage, or penetration of any spinal bone (including the sacrum, pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site at, above, and/or below the level of surgery.
- Retropulsed graft.
- Herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery.
- Non-union (or pseudarthrosis), delayed union, and mal-union.
- Cessation of any potential growth of the operated portion of the spine.
- Loss of or increase in spinal mobility or function.
- Inability to perform the activities of daily living.
- Bone loss or decrease in bone density, possibly caused by stresses shielding.

Important Information for CD HORIZON® Spinal System *continued*

- Graft donor site complications including pain, fracture, or wound healing problems.
- Ileus, gastritis, bowel obstruction, loss of bowel control, or other types of gastrointestinal system compromise.
- Hemorrhage, hematoma, occlusion, seroma, edema, hypertension, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, damage to blood vessels, or other types of cardiovascular system compromise.
- Reproductive system compromise, including sterility, loss of consortium, and sexual dysfunction.
- Development of respiratory problems, e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.
- Change in mental status.
- Death.

Note: Additional surgery may be necessary to correct some of these potential adverse events.

ADDITIONAL POTENTIAL ADVERSE EVENTS FOR PEDIATRIC PATIENTS

- Inability to use pedicle screw fixation due to anatomic limitations (pedicle dimensions, distorted anatomy)
- Pedicle screw malpositioning, with or without neurological or vascular injury
- Proximal or distal junctional kyphosis
- Pancreatitis

WARNING

The safety and effectiveness of pedicle screw spinal systems have been established only for spinal conditions with significant mechanical instability or deformity requiring fusion with instrumentation. These conditions are significant mechanical instability or deformity of the thoracic, lumbar, and sacral spine secondary to degenerative spondylolisthesis with objective evidence of neurologic impairment, fracture, dislocation, scoliosis, kyphosis, spinal tumor, and failed previous fusion (pseudarthrosis). The safety and effectiveness of this device for any other conditions are unknown. The implants are not prostheses. In the absence of fusion, the instrumentation and/or one or more of its components can be expected to pull out, bend, or fracture as a result of exposure to every day mechanical stresses.

A device that has been implanted should never be reused, reprocessed or resterilized under any circumstances. Sterile packaged devices should also never be resterilized. Reuse, reprocessing, or resterilization may compromise the structural integrity of these implants and create a risk of contamination of the implants which could result in patient injury, illness, or death.

ADDITIONAL WARNINGS FOR PEDIATRIC PATIENTS

Warning: The safety and effectiveness of this device has not been established for use as part of a growing rod construct. This device is only intended to be used when definitive fusion is being performed at all instrumented levels.

The use of pedicle screw fixation in the pediatric population may present additional risks when patients are of smaller stature and skeletally immature. Pediatric patients may have smaller spinal structures (pedicle diameter or length) that may preclude the use of pedicle screws or increase the risk of pedicle screw malpositioning and neurological or vascular injury. Patients not skeletally mature that undergo spinal fusion procedures may have reduced longitudinal spinal growth, or may be at risk for rotational spinal deformities (the "crankshaft phenomenon") due to continued differential growth of the anterior spine.

Other adverse events related to pedicle screw fixation, such as screw or rod bending, breakage, or loosening, may also occur in pediatric patients. Pediatric patients may be at increased risk for device-related injury because of their smaller stature.

ADDITIONAL WARNING FOR THE CD HORIZON® SPIRE™ SPINOUS PROCESS PLATE

Please consider the extent of decompression, as well as the amount of intact bone remaining on the spinous processes, when using the CD HORIZON® SPIRE™ Plate as the sole supplemental fixation for an interbody fusion procedure.

PRECAUTIONS

The implantation of pedicle screw spinal systems should be performed only by experienced spinal surgeons with specific training in the use of this pedicle screw spinal system because this is a technically demanding procedure presenting a risk of serious injury to the patient.

A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where many extenuating circumstances may compromise the results. This device system is not intended to be the sole means of spinal support. Use of this product without a bone graft or in cases that develop into a non-union will not be successful. No spinal implant can withstand body loads without the support of bone. In this event, bending, loosening, disassembly, and/or breakage of the device(s) will eventually occur.

Preoperative and operating procedures, including knowledge of surgical techniques, good reduction, and proper selection and placement of the implants are important considerations in the successful utilization of the system by the surgeon. Further, the proper selection and compliance of the patient will greatly affect the results. Patients who smoke have been shown to have an increased incidence of non-unions. These patients should be advised of this fact and warned of this consequence. Obese, malnourished, and/or alcohol abuse patients are also poor candidates for spine fusion. Patients with poor muscle and bone quality and/or nerve paralysis are also poor candidates for spine fusion.

ADDITIONAL PRECAUTIONS FOR PEDIATRIC PATIENTS

The implantation of pedicle screw spinal systems in pediatric patients should be performed only by experienced spinal surgeons with specific training in the use of this pedicle screw spinal system in pediatric patients because this is a technically demanding procedure presenting a risk of serious injury to the patient.

Preoperative and operating procedures, including knowledge of surgical techniques, good reduction, and proper selection and placement of the implants are important considerations in the successful utilization of the system in pediatric patients.

The selection of the proper size, shape, and design of the implant for each patient is crucial to the safe use of this device in pediatric patients.

PHYSICIAN NOTE: Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.



For US Audiences Only

CAUTION: FEDERAL LAW (USA) RESTRICTS THESE DEVICES TO SALE BY OR ON THE ORDER OF A PHYSICIAN.



Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel: + 31 45 566 80 00



Medtronic Sofamor Danek USA, Inc.
1800 Pyramid Place
Memphis, TN 38132
Telephone 800 933 2635 (In U.S.A.)
901 396 3133 (Outside of U.S.A.)
Fax 901 396 0356

©2014 Medtronic Sofamor Danek USA, Inc. All rights reserved.

Please contact your Sales Representative or Customer Service for the most up-to-date version of the package insert.

Summary of Indications

PROGENIX® Putty and PROGENIX® Plus are intended for use as bone graft substitute for voids or gaps that are not intrinsic to the stability of the bony structure (i.e., spine, pelvis, and extremities). It is intended for treatment of surgically or traumatically created osseous defects. Additionally, PROGENIX® Putty and PROGENIX® Plus are intended for the augmentation of deficient maxillary and mandibular alveolar ridges and the treatment of oral/maxillofacial and dental intraosseous defects. Both PROGENIX® Putty and PROGENIX® Plus can be mixed with autograft. When used in spine, PROGENIX® Putty must be mixed with autograft.

MASTERGRAFT® Ceramic Scaffolds are cleared as bone void fillers for bony voids of the skeletal system (i.e., posterolateral spine, pelvis, ilium, and/or extremities). MASTERGRAFT® Granules, MASTERGRAFT® Putty, and MASTERGRAFT® Strip are also cleared as autogenous bone graft extenders. MASTERGRAFT® Granules, MASTERGRAFT® Mini Granules, and MASTERGRAFT® Putty are also cleared for use in bony voids or gaps to fill and/or augment dental oral/maxillofacial bony tissue.

The CAPSTONE® Spinal System is indicated for interbody fusion with autogenous bone graft in patients with Degenerative Disc Disease (DDD) at one or two levels from L2 to S1. These DDD patients may also have up to Grade 1 Spondylolisthesis or retrolisthesis at the involved levels. DDD is defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies. These patients should be skeletally mature and have had six months of non-operative treatment. These implants may be implanted via an open or a minimally invasive posterior approach. Alternatively, these implants may also be implanted via an anterior and/or transforaminal approach. These implants are to be used with autogenous bone graft. These devices are intended to be used with supplemental fixation instrumentation, which has been cleared by the FDA for use in the lumbar spine.

The CLYDESDALE® Spinal System is designed to be used with autogenous bone graft to facilitate interbody fusion and is intended for use with supplemental fixation systems cleared for use in the lumbar spine. The CLYDESDALE® Spinal System is used for patients diagnosed with Degenerative Disc Disease (DDD) at one or two contiguous levels from L2 to S1. These DDD patients may also have up to Grade 1 Spondylolisthesis or retrolisthesis at the involved levels. DDD is defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies. These patients should be skeletally mature and have had six months of non-operative treatment. These implants may be implanted via a minimally invasive lateral approach.

The SOVEREIGN® Spinal System is indicated for use with autogenous bone graft in patients with degenerative disc disease (DDD) at one or two contiguous levels from L2 to S1. DDD is defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies. These patients should be skeletally mature and have had six months of non-operative treatment. These implants may be implanted via a laparoscopic or an open anterior approach. The SOVEREIGN® interbody system may be used as a stand-alone device or in conjunction with supplemental fixation. When used as a stand-alone device, the SOVEREIGN® interbody device is intended to be used with the three titanium alloy fixed or variable angle screws. The accompanying cover plate **MUST** be used anytime the device is used with any number of variable angle screws. If the physician chooses to use less than three or none of the provided screws, then additional supplemental fixation for use in the lumbar spine must be used to augment stability.

The CRESCENT® Spinal System is indicated for interbody fusion with autogenous bone graft in patients with Degenerative Disc Disease (DDD) at one or two levels from L2 to S1. These DDD patients may also have up to Grade 1 spondylolisthesis or retrolisthesis at the involved levels. DDD is defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies. These patients should be skeletally mature and have had six months of non-operative treatment. These implants are to be used with autogenous bone graft. These devices are intended to be used with supplemental fixation instrumentation which has been cleared by the FDA for use in the lumbar spine.

The T2 ALTITUDE™ Expandable Corpectomy System is a vertebral body replacement system intended for use in the thoracolumbar spine (T1-L5) to replace a collapsed, damaged, or unstable vertebral body due to tumor or trauma (i.e., fracture). The T2 ALTITUDE™ Expandable Centerpiece may be used with or without optional modular endcaps which accommodate individual anatomic requirements. The device is to be used with supplemental fixation. Specifically, the construct is to be used with the VANTAGE Anterior Fixation System, the TSRH® Spinal System, the CD HORIZON® Spinal System, or their successors. Additionally, the T2 ALTITUDE™ Expandable Corpectomy System is intended to be used with allograft and/or autograft.

www.medtronic.com

Medtronic
Spinal and Biologics Business
Worldwide Headquarters

2600 Sofamor Danek Drive
Memphis, TN 38132



Medtronic Sofamor Danek USA, Inc.
1800 Pyramid Place
Memphis, TN 38132

(901) 396-3133
(800) 876-3133
Customer Service: (800) 933-2635

www.medtronic.eu

Medtronic International Trading Sàrl
Case postale
Route du Molliau 31
CH-1131 Tolochenaz

Tel: +41 (0)21 802 70 00
Fax: +41 (0)21 802 79 00

The surgical technique shown is for illustrative purposes only. The technique(s) actually employed in each case will always depend upon the medical judgment of the surgeon exercised before and during surgery as to the best mode of treatment for each patient.

Please see the package insert for the complete list of indications, warnings, precautions, and other important medical information.



Consult instructions for use at this website www.medtronic.com/manuals.

Note: Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser.

