



CoreLink.
The Source for Spine™

INSIGHT | PERFORMANCE | VALUE

F3D

INTERBODY SYSTEMS

Featuring patented Mimetic Metal®
Titanium Alloy 3D printing technology



100% Open-pore Architecture



Hydrophilic Wicking




Designed to Reduce Stress Shielding



THE F3D ADVANTAGE

CoreLink's F3D interbodies leverage the benefits of Mimetic Metal which emulates key characteristics of natural bone to provide an optimal structure and environment for healing.



F3D PRODUCT SPECIFICATIONS

Special Order options included. Contact CoreLink Customer Service for more information.

CERVICAL

- **Footprints:** 14.5mm x 12mm
16.5mm x 14mm
- **Heights:** 5mm–12mm
(1mm increments)
- **Lordosis:** 0°, 7°



CORPECTOMY

- **Footprints:** 14.5mm x 12mm
16.5mm x 14mm
- **Heights:** 14mm–50mm
(2mm increments)
- **Lordosis:** 0°, 7°



STRAIGHT

- **Lengths:** 22mm, 26mm, 30mm, 34mm
- **Heights:** 6mm–18mm
(1mm increments)
- **Lordosis:** 0°, 7°, 12°



ANTERIOR LUMBAR

- **Footprints:** 30mm x 22mm
35mm x 25mm
40mm x 27mm
45mm x 30mm
- **Heights:** 10mm–21mm
(1mm increments)
- **Lordosis:** 8°, 15°




LATERAL

- **Widths:** 18mm, 22mm, 25mm
- **Lengths:** 35mm–60mm
(5mm increments)
- **Heights:** 8mm–16mm
(2mm increments)
- **Lordosis:** 0°, 8°, 12°



CURVED


- **Lengths:** 28mm, 31mm
- **Heights:** 7mm–15mm
(1mm increments)
- **Lordosis:** 0°, 6°, 12°



F3D INTERBODY SYSTEMS FEATURE PATENTED MIMETIC METAL TECHNOLOGY

TRABECULAR STRUCTURE


- Designed to emulate key characteristics of natural bone
- Open and randomized pores provide channels for enhanced blood flow to promote fusion



500 MICRON PORES

DYNAMIC FRAMEWORK

- Enables microdeflection to provide more load sharing
- Greater bony on-growth and in-growth compared to machined titanium shown in an in-vivo sheep model*



HISTOLOGY AT 12 WEEKS

*Data on file, pre-clinical data may not be representative of clinical results.

For product information, including indications for use, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the product Instructions for Use (IFU) and information on the CoreLink website. To obtain a copy of the current Instructions for Use (IFU) for full prescribing and risk information, please call CoreLink Customer Service at (888) 349-7808 or visit corelinksurgical.com/ifu. © 2023 CoreLink, LLC. All Rights Reserved. CoreLink is part of CoreLink Holdings, LLC. U.S. Patent No. 10,512,545.



Designed, developed, and manufactured
in the United States of America.

St. Louis, MO
CoreLinkSurgical.com

P: (888) 349-7808
F: (888) 629-5715

CL-MKTG-0070-Rev. 7