

BIO Wedge

Iliac Crest

Description

The BIO Wedge Iliac Crest allografts are recovered from the crest from the ilium, and designed for use in anterior cervical procedures, corpectomy and ankle fusions.

Technical specifications

- Minimum 25mm anterior-posterior depth
- Minimum 8mm medial-lateral width
- 6-10mm and 18mm heights
- Parallel

Product features¹⁸

- Preservon treated for fully hydrated ambient temperature storage
- Structural support provided by cortical aspects
- Provides an osteoconductive matrix for incorporation
- Pre-sized for convenience and intraoperative efficiency
- Processed with Allowash XG to achieve a Sterility Assurance Level to 10^{-6} (SAL)²⁰



BIO Wedge Iliac Crest

BIO Wedge Iliac Crest

Reference number	Description (height)*
7770806	6H
7770807	7H
7770808	8H
7770809	9H
7770810	10H
7770818	18H

*height measurements in mm