




CoreLink.
The Source for Spine™

INSIGHT | PERFORMANCE | VALUE

TIGER® OCT

SPINE FIXATION SYSTEM



COMPREHENSIVE SYSTEM

With a wide assortment of instruments and implants, our system can be customized for preferred surgical technique and individual comfort.

EXCEPTIONAL INSTRUMENTATION

Our instruments are designed for efficiency and ease of use and are crafted in the U.S. by our expert technicians.

THE TIGER OCT ADVANTAGE

The Tiger OCT System provides a full complement of instrument and implant options that seamlessly integrate from occiput to thoracic spine and tie into the Tiger Posterior Lumbar System. The system consists of occipital plates, occipital bone screws, cervical hooks, cervical and upper thoracic bone screws, titanium and cobalt chrome rods, transverse connectors, transverse connectors, and other assorted connectors.

PRODUCT SPECIFICATIONS

SCREWS

Biased Angle Screws (Medial/Lateral and Cephalad/Caudal)


3.5mm: 10mm – 26mm
4.0mm: 10mm – 26mm
4.5mm: 14mm – 26mm

Polyaxial Screws

3.5mm: 10mm – 28mm
4.0mm: 10mm – 28mm
4.5mm: 18mm – 50mm

Long Shank Screws

3.5mm: 24mm – 40mm
4.0mm: 24mm – 40mm




OCCIPITAL PLATE

- Plate sizes: Small, Medium, Large
- Rotating/translating rod connection allows ideal rod alignment
- Five points of bone fixation



Occipital screws:

4.0mm: 6mm – 16mm (2mm increments)
5.0mm: 6mm – 16mm (2mm increments)



4.0MM




5.0MM

*Additional sizes may be available upon request

HOOKS AND CONNECTORS

Cervical Hooks

- Four cervical hook options to best suit patient anatomy



5MM
SMALL

5MM
MEDIUM



5MM
LARGE

7MM

Transverse Connectors

- Head-to-Head and Rod-to-Rod in various lengths



ROD-TO-ROD
CONNECTOR



RETENTION
NUTS




HEAD-TO-HEAD
CONNECTOR




LONG SET
SCREWS

Rod-to-Rod Connectors

- End-to-End and Side-by-Side options
- 3.5mm to 3.5mm or 5.5mm rod




3.5MM TO 5.5MM
SIDE-BY-SIDE
CLOSED



3.5MM TO 5.5MM
SIDE-BY-SIDE
OPEN

Lateral Connectors

- Top Loading and Closed options
- Lengths: 8.5mm, 16mm, 30mm



TOP
LOADING

CLOSED

For product information, including indications for use, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the product Instructions For Use (IFU) and information on the CoreLink website. To obtain a copy of the current Instructions for Use (IFU) for full prescribing and risk information, please call CoreLink Customer Service at (888) 349-7808 or visit corelinksurgical.com/ifu.

© 2020 CoreLink, LLC. All rights reserved.



Designed, developed, and manufactured
in the United States of America.

St. Louis, MO 63026
CoreLinkSurgical.com

P: (888) 349-7808

CL-MKTG-0024, Rev. 2