




CoreLink.
The Source for Spine™

INSIGHT | PERFORMANCE | VALUE

CENTRAFIX™

MIDLINE FIXATION SYSTEM



INCREASED VISUALIZATION

Modular tulip heads allow for intraoperative flexibility and increased visualization with smaller incisions

STRENGTH WITHOUT COMPROMISE

Ultra-low-profile tulip heads accept 4.75mm and 5.5mm rods


THE CENTRAFIX ADVANTAGE

The CoreLink CentraFix Midline Fixation System is a posterior thoracolumbar pedicle screw system designed for a midline approach with a medial to lateral cortical bone trajectory. The system contains modular tulip heads and screw shanks designed specifically for cortical fixation. CentraFix allows for unmatched intraoperative visualization and surgical flexibility.

PRODUCT SPECIFICATIONS


Tulip Heads and Rods:

- Cobalt-Chrome-Molybdenum (CoCrMo) alloy tulip heads accept 4.75mm or 5.5mm diameter rod options
- 4.75mm diameter CoCrMo and 5.5mm diameter* CoCrMo or Ti6Al4V rods
- Straight and pre-lordosed configurations



Screws:

- Diameters: 4.5mm, 5.5mm, 6.5mm, 7.5mm*
- Lengths: 25mm – 45mm (5mm increments)
- Self-drilling-self-tapping tip design



Note: Lengths vary by diameter size. Additional sizes may be available upon request. *Special Order

FEATURES

- Low-profile, anatomic tulip heads allow for reduced tissue disruption without compromising strength
- Modular tulip heads allow for increased visibility and simplified distraction
- Self-drilling, self-tapping cortical screw threading provides easy screw starting in the intended trajectory
- Set screws designed to minimize tulip splay and reduce cross-threading
- Screws allow for 360° motion with 60° cone of angulation



For product information, including indications for use, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the product Instructions For Use (IFU) and information on the CoreLink website. To obtain a copy of the current Instructions for Use (IFU) for full prescribing and risk information, please call CoreLink Customer Service at (888) 349-7808 or visit corelinksurgical.com/ifu.

© 2021 CoreLink, LLC. All rights reserved.



Designed, developed, and manufactured
in the United States of America.

St. Louis, MO 63026
CoreLinkSurgical.com

P: (888) 349-7808

CL-MKTG-0149, Rev. 1