


PRODIGY[®]

MIS PEDICLE SCREW SYSTEM

Features & Benefits




CURITEVA[®]

PRODIGY[®] MIS Pedicle Screw System

The Prodigy MIS Pedicle Screw System is designed to minimize the complexity of procedural steps while correcting and stabilizing the posterior thoracolumbar spine through a percutaneous approach. The system is composed of pedicle screws in standard and dual-thread cortico-cancellous options, ranging from 5.5mm to 7.5mm in diameter in 30 -55mm lengths. The Prodigy MIS system provides reliable rod delivery and significant rod reduction capability while minimizing the incision size and tissue dissection.

MIS Polyaxial Screw

- Color coded tulip heads by bone screw diameter
- Easy tab removal leaving a low profile and clean break
- Simple and efficient instruments provide advanced manipulation techniques



Prodigy MIS Implant Specification

Pedicle Screw Range	Diameter (mm)	Length (5mm Increments)	Color
Standard Screws	5.5	30-55	Green
	6.5	30-55	Dark Blue
Dual-lead Cortico-Cancellous Screws	7.5	30-55	Teal

MIS Titanium Ø5.5 mm Precontoured Rods	30 - 70 (5mm increments) 70 - 110 (10mm increments)
--	--



For more information or to place an order call: 877.9CURITEVA,
email: customersupport@curiteva.com or visit www.curiteva.com

Products patented and/or patent pending.
All products are not currently available in all markets.

Copyright©2021 by Curiteva, Inc.
QMSF-1179 Rev 1.0