

A COMPREHENSIVE SYSTEM FOR TRANSFORAMINAL LUMBAR INTERBODY FUSION



KEY FEATURES & BENEFITS

MULTIPLE FOOTPRINTS

- Three lordotic footprint options conform to various patient anatomies

ANATOMIC DESIGN

- Implant mirrors the natural shape of the intervertebral disc anatomy

TAPERED LEADING EDGE

- Facilitates ease of insertion

AGGRESSIVE TOOTH PATTERN

- Provides stability and reduces the risk of implant migration

GRAFT FENESTRATIONS

- Allow for a copious amount of bone grafting material which may result in increased bony ingrowth

MATERIAL PROPERTIES

- Cages are manufactured from PEEK-OPTIMA®, a biocompatible, radiolucent polymer, which enhances visualization of bone growth in and around the implant
- Embedded tantalum markers facilitate intraoperative placement and radiographic confirmation of implant position

OPERATIVE FLEXIBILITY

- Articulating inserter allows for precise positioning of implant



ANYPLUS® TLIF IMPLANT SIZING

A versatile system that enables a bilateral, unilateral or transforaminal approach to the vertebral bodies.

Simple and intuitive instrumentation facilitates cage insertion.

		9mm Width										HEIGHT (mm)	DEGREE
		7	8	9	10	11	12	13	14	15	16		
LENGTH x WIDTH (mm)	28x9												0°
	31x9												0°
	28x9												4°
	31x9												4°
	28x9												7°
	31x9												7°

		11mm Width										HEIGHT (mm)	DEGREE
		7	8	9	10	11	12	13	14	15	16		
LENGTH x WIDTH (mm)	28x11												0°
	31x11												0°
	28x11												4°
	31x11												4°
	28x11												7°
	31x11												7°

Standard Set

